

LR Reaction

Introduction

An LR reaction inserts one or more parts in pENTR vectors into a pDEST vector. Used to assemble transcriptional units from promoters and genes.

Materials

- ▶ Promoter pENTR plasmid: L4-Promoter-R1
 - ▶ Working concentration: 5 fmol/ul
- ▶ Gene pENTR plasmid: L1-Gene-L2
 - ▶ Working concentration: 5 fmol/ul
- ▶ Destination plasmid: pDEST
 - ▶ Working concentration: 10 fmol/ul
- ▶ Nuclease-free TE
- ▶ 200 μ l PCR strip tubes, 1 tube per rxn
- ▶ 5x LR Clonase II
 - ▶ Stored in ~5 μ l aliquots in the -80 in room 235. **Don't remove an aliquot until you're ready to use it.**
- ▶ Proteinase K
 - ▶ Stored in ~5 μ l aliquots in the -80 in room 235. **Don't remove an aliquot until you're ready to use it.**

Procedure

LR Reaction Setup

1. For each LR you are doing, fill out a column in the following table:

Table1						
	A	B	C	D	E	F
1	Tube Label					
2	Promoter pENTR					
3	Gene pENTR					
4	pDEST					
5						
6						

2. For each LR, label a 200 μ l strip tube with your initials and tube number.

3. Into each tube, pipette:

- 1 μ l of the promoter pENTR
- 1 μ l of the gene pENTR
- 1 μ l of the pDEST

4. Add 1 μ l of TE to each tube

5. Retrieve an aliquot of LR Clonase from the -80.

Bring a razor blade with you, you'll need to cut a tube from the strip tubes.

6. Pulse the LR clonase tube in the microfuge to collect the clonase at the bottom.

7. Add 1 μ l of the LR clonase to each LR reaction.

Be careful pipetting; LR clonase is viscous.

8. Cap the tubes.

9. Flick them several times to mix.

10. Pulse-spin the tubes in the microfuge to collect the liquid at the bottom.

11. Incubate at room temperature for at least 12 hours and not more than 24 hours.

A popular strategy is to tape the tubes to the shelves over the bench, with your initials and the date.

16-24 hours later: Proteinase K kill

12. Retrieve a 5 μ l aliquot of proteinase K from the -80 freezer.

13. Thaw in your fingers, then pulse in the microfuge to collect at the bottom of the tube.

14. Pipette 1 μ l into each of the LR reactions.

15. Flick several times to mix.

16. Pulse-spin the tubes in the microcentrifuge.

17. Incubate at 37° for 15 minutes, or room-temperature for an hour.

PAUSE POINT: You can store the reactions in the -20 indefinitely until the transformation.

00:15:00



18. Proceed to transformation. Transform 2 μ l.

19. Afterwards, cap the tubes. Write the date on the caps and store in the -20 (in case your transformation failed.)

