

Animal Experiments

Construction of hyperuricemia animal models.

- 1) Divide 10 male 5-week-aged balb/c mice into two groups randomly: Group 1 (n=5) and Group 2 (n=5). Set one week's time without special handling to let them adapt to the new environment.
- 2) Feed Group 2 with special feed containing 0.1% adenine since week 2 while feeding Group 1 with normal feed. All the feeding condition is the same expect for feed.
- 3) Take the blood of all the mice from caudal vein after 3 weeks' feeding, isolate the serum and demonstrate the concentration of blood uric acid with ELISA.

Colonization ability test.

- 1) Divide 10 male 5-week-aged balb/c mice into two groups randomly: Group A (n=5) and Group B (n=5). Set one week's time without special handling to let them adapt to the new environment. On the first day of week 2, 200 μ l engineered *E. coli* Nissle 1917 (10^{10} cfu) to Group B, i.g. None special treat is applied to Group A.
- 2) Culture the mice's stools on day 3, 5, 7 and 9 with chloramphenicol which is 3 times of its normal working concentration.
- 3) Detailed steps of stool culturing: take fresh mice's stool into sterilized EP tube, add 0.2ml sterilized ddH₂O, crush it to make turbid liquid. Add 0.8ml sterilized ddH₂O, mix them fully and let stand until the stools sediment. Take 0.1ml supernatant, add LB with chloramphenicol which is 3 times of its normal working concentration, overnight at 37°C.

in vivo test

- 1) Divide 25 male 5-week-aged balb/c mice into five groups randomly: Group Control (n=5), Group HUA (n=5), Group Allopurinol (n=5), Group engineered *E. coli* Nissle 1917 (n=5) and Group *E. coli* Nissle 1917 (n=5). Set one week's time without special handling to let them adapt to the new environment.
- 2) Construct hyperuricemia animal models with all the groups except Group Control.

Apply intragastrical administration as table 1 from week 2, last for 3 weeks.

Group	Feed	i.g. (per day)
Group Control	Normal	200μL ddH ₂ O
Group HUA	Special	200μL ddH ₂ O
Group Allopurinol		200μL allopurinol solution
Group engineered <i>E. coli</i> Nissle 1917		200μL engineered <i>E. coli</i> Nissle 1917 (10 ¹⁰ cfu) on the first day of every week, 200μL ddH ₂ O on other days.
Group <i>E. coli</i> Nissle 1917		200μL <i>E. coli</i> Nissle 1917 (10 ¹⁰ cfu) on the first day of every week, 200μL ddH ₂ O on other days.

Table 1