

## First strand cDNA synthesis

Adapted from: <https://www.neb.com/protocols/2012/09/12/first-strand-cdna-synthesis-e5310>

## Aim of the experiment

This experiment can be used to synthesize cDNA with an RNA template using specific primers.

## Materials

- OneTaq RT-PCR kit (NEB E5310S)
- Gene-specific primer (IDT)
- template RNA

## Procedure

1. To a PCR tube add following reagents:

Table 1: Reaction mix

Concentration	Chemicals
up to 2 $\mu$ g	Template RNA
1 $\mu$ M	Primer
fill up to 8 $\mu$ l	nuclease free H <sub>2</sub> O

For every sample, pipette a duplicate (same pipetting scheme) as negative control.

2. Incubate, 5 min, 70 °C to denature RNA.
3. Add 10  $\mu$ l 2x M-MuLV Reaction Mix to both tubes.
4. To one tube, add 2  $\mu$ l M-MuLV Enzyme Mix.

5. Add 2  $\mu\text{l}$  nuclease free  $\text{H}_2\text{O}$  to the other tube. This is the negative control.
6. Incubate, 1 h, 42 °C.
7. Inactivate the enzyme by incubating, 4 min, 80 °C.
8. Store at -20 °C.