

**iGEM TU/e 2017**  
Biomedical Engineering

Eindhoven University of Technology  
Den Dolech 2, 5612 AZ Eindhoven  
The Netherlands  
[2017.igem.org/Team:TU-Eindhoven](http://2017.igem.org/Team:TU-Eindhoven)

## Transformation into NovaBlue

## Table of contents

Transformation into NovaBlue	<b>1</b>	<b>Transformation into NovaBlue</b>	<b>3</b>
	<b>1.1</b>	<b>Materials</b>	<b>3</b>
	<b>1.2</b>	<b>Setup &amp; Protocol</b>	<b>3</b>
	<b>2</b>	<b>Plating of the cells on agar plates</b>	<b>3</b>

# 1 Transformation into NovaBlue

**Estimated bench time:** 30 minutes

**Estimated total time:** 1.5 hour

**Purpose:** Amplification of the ligation product

It is essential to work sterile, thus disinfect your hands and work near a Bunsen Burner.

## 1.1 Materials

- Bucket with ice
- Bunsen Burner
- Eppendorf tubes with 20 ml competent NovaBlue cells
- Heat/shaking-block
- Incubator
- LB-agar plates supplemented with the correct antibiotic
- Ligation mixture
- Pipettes and tips
- SOC solution (Super optimal broth with catabolite repression)
- Water bath

## 1.2 Setup & Protocol

- Switch on the water bath and set temperature at 42 °C. Also turn on the heat/shaking-block and set up to 37 °C.
- Load a bucket with ice from the ice machine.
- Take the bacterial cells and SOC out of the -80 °C freezer. Transfer the cells directly to ice. Do not touch the bottom of the tube that contains the cells.
- Thaw the cells on ice for ~5 minutes.
- Add 1 µl of ligation mixture to 20 µl bacteria. Leave the bacteria on ice. Mix well. Make sure you work near the Bunsen burner flame.
- Incubate on ice for 5 min.
- Heat shock the cells for exactly 30 seconds at 42 °C.
- Return the cells directly to ice for 2 minutes.
- Add 80 µl of SOC solution (room temperature) to the bacteria. Do not return to ice!
- Incubate for variable time at 37 °C and 300 rpm<sup>1</sup>.
- Let excess water in the LB-agar plates with the correct antibiotics evaporate in the incubator (37 °C) . Place the plates upside down (with the agar up) and slightly opened.

# 2 Plating of the cells on agar plates

**Estimated bench time:** 15 minutes

**Estimated total time:** 18 hour

**Purpose:** Amplification of the ligation product.

<sup>1</sup> Depends on what kind of antibiotic resistance you are working with, see manufacturer protocol for recommended incubation times.

For more information, see our general Plating protocol.