

# 3A Assembly (Digestion and Ligation)

## Aim

To assemble two G-blocks into a new vector.

## Procedure

### Digestion

#### 1. Prepare the following Master Mixes

*Table 1: Master Mixes*

Enzyme Master Mix for Plasmid Backbone (20µl total, for 5 runs)	Enzyme Master Mix for Gblock 1 (BioBrick on the 5' end)	Enzyme Master Mix for Gblock 2 (BioBrick on the 3' end)
2 µl Tango Buffer 1.1 µl EcoRI 1.1 µl PstI 16 µl dH <sub>2</sub> O	2 µl Tango buffer 1.1 µl EcoRI 1.1 µl Kasi 16 µl dH <sub>2</sub> O	2 µl Tango buffer 1.1 µl Kasi 1.1 µl PstI 16 µl dH <sub>2</sub> O

#### 2. Digest Plasmid Backbone

- Add 4 µl linearized plasmid backbone (25 ng/µl for 100 ng total)
- Add 4 µl of Enzyme Master Mix

#### 3. Digest G-block 1

- Add 5 µl G-block 1 (20 ng/µl for 100 ng total)
- Add 5 µl of Enzyme Master Mix

#### 4. Digest G-block 2

**Lab protocol**

Updated: October 28th 2017  
iGEM Stockholm

- Add 5  $\mu$ l G-block 2 (20 ng/ $\mu$ l for 100 ng total)
- Add 5  $\mu$ l of Enzyme Master Mix

5. Digest all three reactions at 37°C for 60 min, heat kill at 80°C for 20 min

## Purification based on the protocol provided by QIAGEN ( QIAquick PCR purification kit)

1. Add 5 volumes Buffer PB to 1 volume of the PCR reaction and mix. (Tips: color  $\rightarrow$  yellow)
2. Place a QIAquick column in a provided 2ml collection tube
3. To bind DNA, apply the sample to the QIAquick column and centrifuge for 30-60 s. Repeat once or twice.
4. To wash, add 600  $\mu$ L buffer PE to the QIAquick column and centrifuge for 30-60 s. Discard flow-through and place the QIAquick column back in the same tube.
5. Centrifuge the QIAquick column once more in the provided 2 ml collection tube for 1 min to remove residual wash buffer.
6. Place each QIAquick column in a clean 1.5 ml microcentrifuge tube
7. To elute DNA, add 50  $\mu$ l buffer EB or water to the center of the QIAquick membrane and centrifuge the column for 1 min. For increased DNA concentration, add 30  $\mu$ l elution buffer to the center of the QIAquick membrane, let the column stand for 1-5 min, centrifuge.
8. If the purified DNA is to be analyzed on a gel add 1 to 5. Mix and Load
9. Measure DNA concentration using Nanodrop

**Lab protocol**

*Updated: October 28th 2017*  
iGEM Stockholm

## Ligation

*Concentration and volumes to be decided on the day.*

1. Add 2µl of digested Plasmid Backbone (50 ng)
2. Add equimolar amount of G-block 1 (EcoRI, KsaI digested) fragment (< 3 µl)\*
3. Add equimolar amount of G-block 2 (KsaI, PstI digested fragment) (< 3 µl)\*
4. Add 1 µl T4 DNA ligase buffer
5. Add 0.5 µl T4 DNA ligase
6. Add water to 10 µl
7. Ligate at RT for 10 min, heat kill at 80°C for 20 min
8. Transform with 1-2 µl of product

\*Easiest to calculate using

[http://www.insilico.uni-duesseldorf.de/Lig\\_Input.html](http://www.insilico.uni-duesseldorf.de/Lig_Input.html)

Start by using a 1:3 vector:insert ratio. If there is enough digested backbone and G-Block left, it is advised to try additional ratios.

## NOTE!

*We used exceedingly more DNA than iGEM did - DNA volumes of 500-1000  $\mu$ L.*

### Recommended Protocol for digestion for KasI (ThermoFisher)

#### 1. Add:

*Nuclease-free water 16  $\mu$ L*

*10X Buffer Tango 2  $\mu$ L*

*DNA (0.5-1  $\mu$ g/ $\mu$ L) 1  $\mu$ L and*

*SspDI 0.5-2  $\mu$ L*

2. Mix gently and spin down for a few seconds.

3. Incubate at 37°C for 1-16 hours.

The digestion reaction may be scaled either up or down.

## Sources

This protocol is modified from the 3A Assembly protocol at [igem.org](http://parts.igem.org/Help:Assembly/3A_Assembly):

[http://parts.igem.org/Help:Assembly/3A\\_Assembly](http://parts.igem.org/Help:Assembly/3A_Assembly)

ThermoFisher

<https://www.thermofisher.com/order/catalog/product/ER2191>

**Lab protocol**

*Updated: October 28th 2017*

**iGEM Stockholm**