

July 31, 2017


Title: Testing for Optimal Fluoride Levels and Dilutions for CHOP Growth

Expected Outcome: The increase of fluoride will increase the growth of CHOP, until a certain point. This experiment will ascertain that point.

Procedure:

- I. Make 1 set of plates with no fluoride, 1 μ M fluoride, 25 μ M fluoride, 50 μ M fluoride, 75 μ M fluoride, 100 μ M fluoride, 150 μ M fluoride. Also add 50 μ g/mL chloramphenicol.
 - A. Make serial dilutions of CHOP and Δ CRCB, and plate them
 - B. Incubate the plates at 37°C and image daily
- II. Make 1 set of plates with no fluoride, 1 μ M fluoride, 25 μ M fluoride, 50 μ M fluoride, 75 μ M fluoride, 100 μ M fluoride, 150 μ M fluoride. Also add 50 μ g/mL chloramphenicol.
 - A. Divide the plates into 4 quadrants.
 - B. Make CHOP at a 10^{-5} dilution and 10^{-4} dilution
 - C. On the bottom half of the plates, use a glass spreader to apply the 10^{-5} dilution
 - D. Apply the 10^{-4} dilution on the top half using a glass spreader
 - E. Incubate the plates at 37°C and image daily

Results:



Conclusion: Based off the plate images, CHOP is best thrives the most at 100 μ M fluoride and at a 10^{-4} dilution.