

DNA kit plate

Aim

To extract biobrick plasmids from the 2017 iGEM Distribution Kit.

Procedure

Note: There is an estimated 2-3ng of DNA in each well, following this protocol, assume that you are transforming with 200-300pg/ μ L

1. With a pipette tip, punch a hole through the foil cover into the corresponding well of the part that you want. Make sure you have properly oriented the plate. Do not remove the foil cover, as it could lead to cross contamination between the wells.
2. Pipette 10 μ L of dH₂O (distilled water) into the well. Pipette up and down a few times and let sit for 5 minutes to make sure the dried DNA is fully resuspended. The resuspension will be red, as the dried DNA has cresol red dye. We recommend that you do not use TE to resuspend the dried DNA.
3. Transform 1 μ L of the resuspended DNA into your desired competent cells, plate your transformation with the appropriate antibiotic* and grow overnight.
4. Pick a single colony and inoculate broth (again, with the correct antibiotic) and grow for 16 hours.
5. Use the resulting culture to miniprep the DNA AND make your own glycerol stock (for further instruction on making a glycerol see this page). We recommend using the miniprep DNA to run QC tests, such as restriction digests and sequencing.

Lab protocol

Updated: October 28th 2017
iGEM Stockholm

Sources

This protocol is a modified version of a protocol by Gunilla Karlsson Hedestam Group at Karolinska Institutet.

Lab protocol

Updated: October 28th 2017

iGEM Stockholm