

USER linearization

To linearize vectors by opening the USER cassette, we use the following digestions.

The first digestion cuts the cassette, the second creates the overhang that will be used for insertion of the genes in the vector.

Protocol

1. Digestion 1

Component	100 µl reaction
pET102 iGEM	42 µl
AsiSI enzyme	2.5 µl
Thermofisher tango buffer	10 µl
Nuclease free water	45.5 µl

2. Mix and incubate for 3 hours at 37 °C

3. Column purify (gel extraction protocol, but add 150 µl binding buffer no matter the amount of gel)

4. Digestion 2

Component	
Elution from purification	50 µl
Nb.BsmI enzyme	1 µl
NEBuffer 3.1	6 µl
Nuclease free water	3 µl