



MATERIALS

PROTOCOL FOR CHARACTERIZATION OF THE SENSING MODULE

1X PBS

AHL SOLUTION

9 DIFFERENT PARTS LB+CHL

96 DEEP WELL PLATE

96 WELL PLATE

DAY 1

1. Transform the followings into DH10B E.Coli:
 - a. pSB1C3_BBa_F2620 (-)
 - b. pSB1C3_BBa_J23110_E0240 (+)
 - c. pSSB1C3_BBa_J23117_E0240 (+)
 - d. pSB1C3_BBa_T9002
 - e. pSB1C3_BBa_F2620-C0261-E0240
 - f. pSB1C3_BBa_F2620-C0261-E0240-Anti1
 - g. pSB1C3_BBa_F2620-C0261-E0240-Anti2
 - h. pSB1C3_BBa_F2620-ABR1-C0261-ABR1-E0240-Anti1
 - i. pSB1C3_BBa_F2620-ABR2-C0261-ABR2-E0240-Anti2

DAY 2

2. Pick a colony on each plate and streak the following:
 - a. pSB1C3_BBa_F2620 (-)
 - b. pSB1C3_BBa_J23110_E0240 (+)
 - c. pSSB1C3_BBa_J23117_E0240 (+)
 - d. pSB1C3_BBa_T9002
 - e. pSB1C3_BBa_F2620-C0261-E0240
 - f. pSB1C3_BBa_F2620-C0261-E0240-Anti1
 - g. pSB1C3_BBa_F2620-C0261-E0240-Anti2
 - h. pSB1C3_BBa_F2620-ABR1-C0261-ABR1-E0240-Anti1
 - i. pSB1C3_BBa_F2620-ABR2-C0261-ABR2-E0240-Anti2

DAY 3

1. Pick 3 colonies from each streaked plates for inoculation in 5mL LB+CHL (there should be a total of $9 \times 3 = 27$ tubes):
 - a. pSB1C3_BBa_F2620 (-)
 - b. pSB1C3_BBa_J23110_E0240 (+)
 - c. pSSB1C3_BBa_J23117_E0240 (+)
 - d. pSB1C3_BBa_T9002
 - e. pSB1C3_BBa_F2620-C0261-E0240
 - f. pSB1C3_BBa_F2620-C0261-E0240-Anti1
 - g. pSB1C3_BBa_F2620-C0261-E0240-Anti2
 - h. pSB1C3_BBa_F2620-ABR1-C0261-ABR1-E0240-Anti1
 - i. pSB1C3_BBa_F2620-ABR2-C0261-ABR2-E0240-Anti2

DAY 4

1. Sample dilution:
 - a. Measure the OD600 of each culture.
 - b. Dilute the overnight culture so that it becomes 8ml 0.1 OD600 culture.
 - c. Incubate them inside the shaker until it become 0.6 (around 2 hours, but start measuring OD600 after 1hr 30mins).

Meanwhile,

2. AHL Preparation:
 - a. Make 4.38×10^{-4} M AHL stock solution by AHL powder and DMSO.
 - b. Prepare 20 tubes of Eppendorf.
 - c. Transfer 60ul 4.38×10^{-4} M AHL into the first Eppendorf, and make 19 tubes of 54ul 1XPBS.
 - d. Make a serial dilution from 4.38×10^{-4} M to 4.38×10^{-23} M by transfer 6ul previous AHL solution into 54ul 1XPBS. *Make sure in each dilution, pipette up and down several times to mix well.
3. Sample loading after OD600 reaches 0.6:
 - a. For 96 Deep Well plate: (3 plates)
 - i. pSB1C3-BBa_T9002: transfer 11 x 650ul of cell culture into for A1-11
 - ii. pSB1C3-BBa_F2620-C0261-E0240: transfer 11 x 650ul of cell culture into for B1-11
 - iii. pSB1C3-BBa_F2620-C0261-E0240-R0063-Anti1: transfer 11 x 650ul of cell culture into for C1-11
 - iv. pSB1C3-BBa_F2620-C0261-E0240-R0063-Anti2: transfer 11 x 650ul of cell culture into for D1-11
 - v. pSB1C3-BBa_F2620-C0261-E0240(ABR1)-R0063-Anti1: transfer 11 x 650ul of cell culture into for E1-11



- vi. pSB1C3-BBa_F2620-C0261-E0240(ABR2)-R0063-Anti2: transfer 11 x 650ul of cell culture into for D1-11
- vii. pSB1C3-BBa_J23110-E0240: transfer 1 x 650ul of cell culture into for H1
- viii. pSB1C3-BBa_J23117-E0240: transfer 1 x 650ul of cell culture into for H2
- ix. pSB1C3-BBa_F2620: transfer 1 x 650ul of cell culture into for H3
- x. LB+CHL: transfer 1 x 650ul of cell culture into for H4

- b. For each part and each colony, transfer 100ul of cell culture into each well and make triplicate (3x100ul) at 0min. (See Appendix 1) (3x3=9 plates)
- c. Then from part ai) to avii), add 8ul of corresponding AHL in the 350ul culture to make 11 different condition solution (0M, 10-5M,10-7M,....,10-23M).

96 deep wells plate:

	1	2	3	4	5	6	7	8	9	10	11	12
A	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
B	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
C	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
D	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
E	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
F	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻²¹ M	10 ⁻²³ M	0M	
G												
H	110	117	F	LB								

[pSB1C3-BBa_T9002: Blue](#)

[pSB1C3-BBa_F2620-C0261-E0240: Green](#)

[pSB1C3-BBa_F2620-C0261-E0240-R0063-Anti1: Light Blue](#)

[pSB1C3-BBa_F2620-C0261-E0240-R0063-Anti2: Blue](#)

[pSB1C3-BBa_F2620-C0261-E0240\(ABR1\)-R0063-Anti: Purple](#)

[pSB1C3-BBa_F2620-C0261-E0240\(ABR2\)-R0063-Anti2: Light Pink](#)

[pSB1C3-BBa_J23110-E0240: Red](#)

[pSB1C3-BBa_J23117-E0240: Orange](#)

[pSB1C3-BBa_F2620: Dark Red](#)

[LB+CHL: Pink](#)

- c. Incubate the 96 deep wells plate into 37C shaker for 3 hrs.
- d. For each part, again, transfer 100ul of cell culture into each well and make triplicate at 200min. (See Appendix 1) (3x3=9 plates)

*All procedures in step 3a) to 3c) and 3e) should be performed on ice.



4. Measurements:

Program Plate Reader:

- a. Shaking: 15 seconds of orbital shaking (6mm amplitude) followed by 10 seconds waiting before the measurement (to ensure consistency within culture sample)
- b. Measure fluorescence(RFU) and OD595
- c. Calculate the Fluroescence/OD
- d. Provide the data with all Raw data to us.

Appendix 1: (Time: 0, 200 mins)

Plate 1: [pSB1C3-BBa T9002: Blue](#), [pSB1C3-BBa F2620-C0261-E0240: Green](#)

	1	2	3	4	5	6	7	8	9	10	11	12
A	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
B	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
C	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
D	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
E	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
F	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
G												
H	110	110	110	117	117	117	F	F	F	LB	LB	LB

Plate 2: [pSB1C3-BBa F2620-C0261-E0240-R0063-Anti1](#), [pSB1C3-BBa F2620-C0261-E0240-R0063-Anti2](#)

	1	2	3	4	5	6	7	8	9	10	11	12
A	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
B	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
C	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
D	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
E	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
F	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
G												
H												

Plate 3: [pSB1C3-BBa F2620-C0261-E0240\(ABR1\)-R0063-Anti1](#), [pSB1C3-BBa F2620-C0261-E0240\(ABR2\)-R0063-Anti2](#)

	1	2	3	4	5	6	7	8	9	10	11	12
A	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
B	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
C	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
D	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁵ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁷ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻⁹ M	10 ⁻¹¹ M	10 ⁻¹¹ M	10 ⁻¹¹ M
E	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹³ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁵ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁷ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M	10 ⁻¹⁹ M
F	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²¹ M	10 ⁻²³ M	10 ⁻²³ M	10 ⁻²³ M	0M	0M	0M			
G												
H												

-----END-----