

### Prodeoxyviolacein *in vitro* assay

(optimized for our application, referring to "In Vitro Biosynthesis of Violacein from L-Tryptophan by the Enzymes VioA–E from *Chromobacterium violaceum*"; Carl J. Balibar and Christopher T. Walsh; Biochemistry, 2006)

The following yeast samples were harvested with 10 total OD<sub>600</sub> units.

WT1	WT2	A1.1	A1.2	A2.1	A2.2	B1.1	B1.2	B2.1	B2.2	E1.1	E1.2	E2.2	E2.2
-----	-----	------	------	------	------	------	------	------	------	------	------	------	------

After this the purification was implemented with a protocol by invitrogen ("pGAPZ A, B, and C pGAPZ  $\alpha$  A, B, and C Pichia expression vectors for constitutive expression and purification of recombinant proteins"; invitrogen by life technologies; user manual, page 23-24, 2010).

This led to two protein- and two cell suspensions of each sample. The replicates were mixed as followed:

WT1 CS	WT2 PS	A1.1 CS	A2.1 CS	A1.2 PS	A2.2 PS	B1.1 CS	B2.1 CS	B1.2 PS	B2.2 PS	E1.1 CS	E2.1 CS	E1.2 PS	E2.2 PS
<b>WT CS</b>	<b>WT PS</b>	<b>VioA CS</b>		<b>VioA PS</b>		<b>VioB CS</b>		<b>VioB PS</b>		<b>VioE CS</b>		<b>VioE PS</b>	

PS = protein suspension, CS = cell suspension

#### mastermix (600 $\mu$ l)

1200 mM glycine (54 mg)

7,8 mM FAD (38 mg)

16 units catalase per 1  $\mu$ l (4,8 mg)

360  $\mu$ l MgCl<sub>2</sub>: 150 mM

240  $\mu$ l ddH<sub>2</sub>O

#### suspension mixtures

WT CS: 1500  $\mu$ l WTZ

WT PS: 1500  $\mu$ l WTP

CS1: 500  $\mu$ l VioA CS, 500  $\mu$ l VioB CS, 500  $\mu$ l VioE CS

PS1: 250  $\mu$ l VioA PS, 1000  $\mu$ l VioB PS, 250  $\mu$ l VioE PS

CS2: 500  $\mu$ l VioA CS, 500  $\mu$ l VioB CS, 500  $\mu$ l VioE CS

PS2: 250  $\mu$ l VioA PS, 1000  $\mu$ l VioB PS, 250  $\mu$ l VioE PS

#### protocol

- each reaction mix: 100  $\mu$ l mastermix, 1468  $\mu$ l suspension mixture
- adjust pH=9,25
- add 32  $\mu$ l 1000  $\mu$ M L-tryptophan to start the reaction
- incubate at 30 degrees; 400 rpm
- take 300  $\mu$ l of each sample at 0, 30, 60, 90, 120 min
- add 600  $\mu$ l icecold MeOH, 60  $\mu$ l icecold DMSO to the sample and mix roughly
- centrifuge 5 min at 13000 rpm
- transfer supernatant in new loBind reaction tube
- incubate 10 min on ice
- centrifuge 5 min at 13000 rpm
- transfer 200  $\mu$ l supernatant into LC-Vial for mass spectrometry analysis