## Biotinylation of glass slides

- Immerse the plates in acetone to remove surface contaminants
- rinse with deionized water and dry under a stream of nitrogen
- ◆ Place in a piranha solution (conc. H₂SO₄/H₂O₂ (2.5:1 v/v)) for 30 min at room temperature
- Wash with deionized water and dry under a stream of nitrogen
- immerse in a freshly prepared solution of APTES (Sigma Aldrich, Australia) and anhydrous toluene (2% (v/v) APTES) for 30 min at room temperature
- rinse with anhydrous toluene to remove any excess reagent
- heat in an oven at 110°C for 20 min
- Fresh solutions of biotin-NHS in DI water with a concentration of 4.0 mg/mL were transferred to the slides
- Seal the slides in a humidity chamber to prevent evaporation of the biotin-NHS solution for 60 min.
- ◆ The slides was then rinsed copiously with DI water and subsequently sonicated in DI water for ~10-15 min.

## References

- Acres, R. G., Ellis, A. V., Alvino, J., Lenahan, C. E., Khodakov, D. A., Metha, G. F., & Andersson, G. G. (2012). Molecular structure of 3-aminopropyltriethoxysilane layers formed on silanol-terminated silicon surfaces. *Journal of Physical Chemistry C*, 116(10), 6289–6297. https://doi.org/10.1021/jp212056s
- Lapin, N. A., & Chabal, Y. J. (2009). Infrared characterization of biotinylated silicon oxide surfaces, surface stability, and specific attachment of streptavidin. *Journal of Physical Chemistry B*, 113(25), 8776–8783. https://doi.org/10.1021/jp809096m

