

iGEM 2017 – Microbiology – BMB – SDU	
Title: Bacterial freezing stock	Date issued: 2012.10.26
SOP number: SOP11	Review date: 2017.05.01
Version number: 02	Written by: MM
	Modified by: JB

1. Purpose

To freeze new strains of bacteria in -80°C

2. Area of application

This procedure is valid for all *E. coli* strains

3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Freezer (-80 °C)	Laboratory (class 1) – V18-404b-0 Laboratory (class 2) – V15-501b-2		-75 °C - -80 °C
Pipette p1000			
15 ml tube rack	Hallway storage (1. Floor)		
Vortex			

4. Materials and reagents – their shelf life and risk labelling

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Glycerol	50 %	AppliChem	Anne Mette, RT	
Blue pipette tips		Contact Lab-manager	Micro storage	
Freezing tubes		Contact Lab-manager	Micro storage	

5. QC – Quality Control

6. List of other SOPs relevant to this SOP

iGEM2017_SOP02_v02_FN_ONC_E.coli

7. Environmental conditions required

8. Procedure

- 8.1 Add 300 µL 50 % glycerol to each sterile freezing tube
- 8.2 Vortex ON culture of bacteria
- 8.3 Add 900 µL ON culture to each freezing tube with glycerol
- 8.4 Pipette up and down to mix
- 8.5 Check register sheet for number for the cultures and write on the tubes
- 8.6 Add to strain collection

9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
ON Culture		Liquid bacterial waste	
Once use plastic		GMO yellow waste	

10. Time consumption

- Total-time 15 min.
- Hands-on-time 15 min.

11. Scheme of development

Date / Initials	Version No.	Description of changes
12.10.25 / MM	01	The SOP has been written
13.01.03 / MM & TK	01	The SOP has been approved
17.05.01 / JB	02	The SOP has been modified
17.05.01 / EG	02	The SOP has been approved

12. Appendices