

# Self-Circularization of Linear DNA

---

## Introduction

This protocol is used to re-circularize linear plasmid DNA. Obtained from UNC Chapel Hill.

## Materials

### › Material

- › 5 uL of 10x T4 DNA ligase buffer
- › 10-50 ng of Linearized Plasmid DNA
- › 5 uL T4 DNA Ligase
- › \*NF water to final volume 50 uL

## Procedure

### Procedure

1. Mix components above -- Start with water, end with Ligase, the order of the other components does not matter
2. Vortex quickly to mix thoroughly
3. Incubate at RT for 2 hours
4. Heat inactivate Ligase at 65 C for 10 mins in thermocycler
5. Transform up to 5 uL of re-circularized DNA into competent cells using the transformation protocol