

Ligation

Aim

Ligation of one or more DNA fragments previously digested with restriction enzymes for further experiments.

Procedure

Initial notes:

The T4 DNA Ligase Buffer should be thawed and resuspended at room temperature.

The T4 DNA Ligase should be added last.

Mix the following components on ice.

| COMPONENT | 20 µl REACTION |
|----------------------------|----------------|
| T4 DNA Ligase Buffer (10X) | 2 µl |
| Vector DNA* | variable |
| Insert DNA* | variable |
| T4 DNA ligase | 1 µl |
| RNase-free water | To 50 µl |

*Molar ratios can be calculated with the [University of Dusseldorf Ligation Calculator](#).

Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube. Follow with a quick ("touch") spin-down in a microcentrifuge. Do not vortex the reaction.

Incubate at 16°C overnight or at room temperature for 10 minutes.

Heat-inactivate at 65°C for 10 minutes.

Lab protocol

Updated: October 28th 2017

iGEM Stockholm

Sources

This protocol is a modified version of the original [Ligation Protocol with T4 DNA Ligase](#) provided by NEB®.

Lab protocol

Updated: October 28th 2017

iGEM Stockholm