

iGEM TU/e 2017
Biomedical Engineering

Eindhoven University of Technology
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
2017.igem.org/Team:TU-Eindhoven

LIU PCR

Table of contents

LUI PCR	1	LIU PCR	3
	1.1	Materials	3
	1.2	Setup & Protocol	3

1 LIU PCR

Estimated bench time: 30 minutes

Estimated total time: 5 hours

Purpose: insertion of amino acids in the DNA sequence with the primers.

It is essential to work with gloves at all times to protect the DNA from DNase activity.

1.1 Materials

- Autoclaved H₂O
- Autoclaved PCR tubes
- Bucket with ice
- DNA polymerase
- dNTPs
- Forward primer
- Pipettes and tips
- Reverse primer
- Thermal cycler
- Vector DNA in which amino acids need to be inserted

1.2 Setup & Protocol

- Construct a PCR mixture in the following way. Start with the component with the largest volume and end with the Master Mix. Keep the Master Mix on ice.

Component	Quantity/mass/final concentration	Volume (µl)
H ₂ O	Fill up to 50 µL	
DNA	110 ng/µL	
Primer FW	1 µM (10 µM stock)	5
Primer RV	1 µM (10 µM stock)	5
Q5 Mastermix (2X)	1X	25
Total		50

- Mix well by pipetting up and down.

Primer	T _m (NO)	T _m (PP)	Concentration
RV_5AA_streptagCT33	52	41	
FW_5AA_streptagCT33	46	41	
Thus for 5AA	44	35	
RV_1AA_streptagCT33	52	39-41	
FW_1AA_streptagCT33	56	39-41	
Thus for 1AA	48	35	

NO=primer sequence matched to the template

PP=Primer-Primer overlap

- Run the following PCR program

Cycle	Temperature	Time	Description
1	95	5 minutes	Denature the template DNA
30x amplification	95	1 minute	
	T_m(NO)-5	1 minute	
	72	10-15 minutes	Depending on the length of the template constructs
PCR finishing cycles	T_m(PP)-5	1 minute	Annealing step
	72	30 minutes	Extension step

- Digestion of PCR product with DpnI. See the protocol Digestion for more information.