

# Protocol for Western blotting

## 1. Preparation of lysate from 60 mm×15 mm cell culture plate:

- ( 1 ) Aspirate the growth medium from the cells, and wash the cells with PBS.
- ( 2 ) Aspirate the PBS, then add 600 $\mu$ L trypsin EDTA solution.
- ( 3 ) Incubation for several minutes at 37 °C and check the culture with inverted microscope to be sure that the cells are rounded up which indicates they are detached from the surface.
- ( 4 ) Add 2 mL fresh complete medium, then pipet all the adherent cells into cell suspension and pipet up and down.
- ( 5 ) Transfer the cell suspension into a pre-cooled 15 mL centrifuge tube and centrifuge at 1200 rpm for 5 minutes.
- ( 6 ) Aspirate the filtrate, add 1mL PBS and resuspend the cell, transfer the cell suspension into a pre-cooled 1.5 mL microcentrifuge tube.
- ( 7 ) Centrifuge at 1200 rpm for 5 minutes at 4 °C, then aspirate the filtrate.
- ( 8 ) Add 35  $\mu$ L RIPA Lysis Buffer containing protease inhibitor cocktail and phospho-stop , and mix or vortex briefly to resuspend the cells completely
- ( 9 ) Incubate on ice for 15 minutes. Vortex briefly to resuspend and lyse residual cells.
- (10) Centrifuge at 13200 rpm for 15 minutes at 4 °C.
- (11) Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice, and discard the pellet.

## 2. Protein Quantification:

We use Beyotime BCA Protein Assay Kit.

- (1) Dilute 10 $\mu$ L BSA Standard (5 mg/mL) to 100 $\mu$ L.
- (2) Prepare Working Reagent by mixing 50 parts of BCA™ Reagent A with 1 part of BCA™ Reagent B (50:1, Reagent A:B).
- (3) a. Pipette 0、 1、 2、 4、 6、 8、 12、 16、 20  $\mu$ L of standard into a 96 microplate well and dilute each standard to 20  $\mu$ L.  
b. Pipette 20  $\mu$ L of each diluted sample replicate into the microplate well.  
c. Add 200  $\mu$ L of the Working Reagent to each well and mix plate thoroughly by pipette.  
d. Cover plate and incubate at 37 °C for 30 minutes.
- ( 4 ) Measure the absorbance at 595 nm on a plate reader.
- ( 5 ) Prepare a standard curve by plotting the average Blank-corrected 595 nm measurement for each BSA standard vs. its concentration in  $\mu$ g/ml. Use the standard curve to determine the protein concentration of each unknown sample.

## 3. Reduce and denature samples:

Boil each cell lysate in LDS sample buffer(4 $\times$ ) at 95 °C for 5 minutes. Lysates can be aliquoted

and stored at -20 °C for future use.

#### 4. Loading and running the gel:

Load 60–80 µg of total protein from cell lysate into the wells of the SDS-PAGE gel(10% Gel percentage), along with pre-stained molecular weight marker.

#### 5. Run the gel for 55 minutes at 180 V.

#### 6. Transfer proteins from the gel to membrane:

(1) a. Prepare transfer buffer(1L): 3.03 g Trizma base ( 25 mM), 14.4 g Glycine ( 192 mM),200 mL methanol(20%),800 ml water;

b. Prepare PBST(1L):PBS(1L), 0.5ml Tween 20 (0.05%)

( 2 ) Cut a piece of PVDF membrane according to the size of gel and wet in methanol.for 2 minutes.

( 3 ) Transfer the membrane to the transfer buffer.

( 4 ) Assemble transfer stack as follows:



( 5 ) Attach the electrodes. Set the power supply to 350 mA (constant current) for 1h 20 min at 4 °C.

#### 7.Blocking:

( 1 ) Prepare Blocking buffer:Add 2.5g non-fat powdered milk ( 5% ) to 50 mL PBST buffer, mix thoroughly.

( 2 ) Remove the membrane to a small container and stain in Ponceau stain reagent for several minutes,then wash for short time with PBST to check the transfer.

( 3 ) Cut the membrane into several stripes according to the size of target protein.

( 4 ) Completely.destain in PBST.

( 5 ) Incubate the membrane in 5% milk for 1 h at room temperature.

#### 8. Primary antibody incubation:

Incubate the membrane with appropriate dilutions of primary antibody in blocking buffer, 4°C

overnight at recommended concentration.

**9. Secondary antibody incubation:**

- ( 1 ) Wash the membrane in 5 washes of PBST, 6 min each.
- ( 2 ) Incubate the membrane with the recommended dilution of conjugated secondary antibody in blocking buffer at room temperature for 1 h.

**10. ECL detection:**

We use Pierce™ ECL Western Blotting Substrate.

- ( 1 ) Wash the membrane in 3 washes of PBST, 5 min each.
- ( 2 ) Prepare the substrate working solution by mixing equal parts of Detection Reagents 1 and 2. Use 0.125mL Working Solution per cm<sup>2</sup> of membrane.
- ( 3 ) Incubate blot with working solution for 1 minute at RT.
- ( 4 ) Remove blot from working solution and place it in a clear plastic wrap. Use an absorbent tissue to remove excess liquid.
- ( 5 ) Place the membrane in a chemiluminescence imaging system with the protein side facing up and carefully press out any bubbles from between the chemiluminescence imaging system and the membrane.
- ( 6 ) Cover the chemiluminescence imaging system and detect.