



## Colony PCR GoTaq® G2 DNA Polymerase

- ◆ First alternative:
  - ◇ Pick one colony with a sterile tip and elute in 100  $\mu\text{L}$  ddH<sub>2</sub>O or medium
  - ◇ Store the colony at 4 °C while colony PCR is running
- ◆ Second alternative:
  - ◇ Pick one colony with a sterile tip and streak cells at marked position on a new plate
  - ◇ Put tip in PCR tube already containing the reaction mixture
- ◆ For both alternatives continue as follows (one reaction mix):
  - ◇ 5  $\mu\text{L}$  5x Green or Colorless GoTaq® Reaction Buffer
  - ◇ 0.5  $\mu\text{L}$  PCR Nucleotide Mix (10 mM each)
  - ◇ 0.25  $\mu\text{L}$  upstream primer (0.1 – 1.0  $\mu\text{M}$ )
  - ◇ 0.25  $\mu\text{L}$  downstream primer (0.1 – 1.0  $\mu\text{M}$ )
  - ◇ 0.25  $\mu\text{L}$  GoTaq® G2 DNA Polymerase (Promega)
  - ◇ 1  $\mu\text{L}$  template DNA
  - ◇ 17.75  $\mu\text{L}$  ddH<sub>2</sub>O

- ◆ PCR program:

	Steps	°C	m:s
	Cell lysis and denaturation	95.0	3:00
30 – 35 x	Denaturation	95.0	0:30
	Primer annealing	annealing temperature	0:30
	Extension	72.0	1 min/kb
	Final extension	72.0	5:00
	Hold	8.0	PAUSE

From: [Promega](#)