

## Native DNA PAGE

- ◆ For each gel (12%) use:
  - ◇ 3.35 mL Bisacrylamid (30 %)
  - ◇ 6.15 mL nuclease-free H<sub>2</sub>O
  - ◇ 0.5 mL 10x TBE buffer
  - ◇ 200 µL APS (10 %)
  - ◇ 10 µL TEMED (99%)
- ◆ Pour the solution quickly into the gel casting form.
- ◆ Insert comb without getting bubbles stuck underneath.
- ◆ Leave the gel at room temperature for >60 minutes to polymerize.
- ◆ Preparing the sample:
  - ◇ 10 µL DNA sample
  - ◇ 4 µL Loading Dye (6x, SDS free)
  - ◇ 6 µL nuclease-free H<sub>2</sub>O
- ◆ Running the gel:
  - ◇ Load the samples and DNA ladder into the gel pockets.
  - ◇ Connect the power lead and run the gel with 60 V for ~2 h (for small fragments, ~80 bp).
- ◆ After finishing the PAGE, remove gel from gel casting form and transfer it into an ethidium bromide bath.