

PCR Using Q5® High-Fidelity 2X Master Mix

Aim

Amplification of a DNA fragment of interest for further experiments.

Procedure

Thaw Q5 High-Fidelity 2X Master Mix, primers, and DNA sample on ice. Let thaw completely and mix before adding to reaction mix to avoid localized concentration differences.

Prepare a primer working solution containing 1 μ l of the forward primer (100 μ M) and 1 μ l of the reverse primer (100 μ M), and 8 μ l of RNase-free water for a final volume of 10 μ l. This achieves a final primer concentration of 200 nM.

Set up the following reaction on ice:

COMPONENT	25 μ l REACTION	50 μ l REACTION	FINAL CONCENTRATION
Q5 High-Fidelity 2X Master Mix	12.5 μ l	25 μ l	1X
Primer mix	0.6 μ l	1.2 μ l	200 nM
Template DNA	variable	variable	<1000 ng
RNase-free water	To 25 μ l	To 50 μ l	

Notes: Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary.

Quickly transfer PCR tubes to a PCR machine and begin thermocycling.

Lab protocol

Updated: October 28th 2017

iGEM Stockholm

STEP	TEMP	TIME
Initial Denaturation	98°C	30s
25-35 cycles: Denaturation Annealing Extension	98°C *50-72°C 72°C	5-10s 10-30s 20-30s/kb
Final Extension	72°C	2-5min
Hold	4°C	-

*Use of the [NEB Tm Calculator](#) is highly recommended.

When the program is finished, the PCR product may be removed and stored at 4°C.

Note:

The above protocol could also be used for colony PCR reactions.

The colony PCR reaction set-up follows the standard PCR protocol described above.

The DNA is obtained by picking a single colony from the plate using a sterile pipette tip and thoroughly swirling it in 20 µl of RNase-free water.

1 µl of this water will be used for the PCR reaction.

The rest could be plated or inoculated in a flask overnight.

Sources

This protocol is a modified version of the original [Q5® High-Fidelity 2X Master Mix PCR amplification protocol](#) provided by NEB®.

Lab protocol

Updated: October 28th 2017

iGEM Stockholm