

iGEM TU/e 2017
Biomedical Engineering

Eindhoven University of Technology
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
2017.igem.org/Team:TU-Eindhoven

Flow chamber

Table of contents

Flow chamber	1	Preparation and loading of the Flow chamber	3
	1.1	Materials	3
	1.2	Setup	3
	1.3	Loading flow chamber	3
	1.3.1	Materials	3
	1.3.2	Setup & Protocol	3

1 Preparation and loading of the Flow chamber

Estimated bench time: 0.5 hours

Estimated total time: depending on the amount of samples roughly 20 mins per sample

Purpose: Preparing samples for microscopy assays.

1.1 Materials

- Cover plate
- BSA (30 mg/ml)
- Double-sided adhesive tape
- Ethanol 100%
- Microscopy glass slides
- samples
- Nail polish
- Tissues fiber-free
- Assay Buffer

1.2 Setup

- Clean the glass plates with 100% ethanol
- Cut the double-sided adhesive tape in half and place two strips of approximately 2-3 cm on the long sides of the glass plates leaving a channel in the middle.
- Clean the cover plate with 100% ethanol and place it over the channel.
- This creates a flow channel with an approx. volume of 20-25 μ l

1.3 Loading flow chamber

1.3.1 Materials

- BSA (30 mg/ml)
- samples
- Nail polish
- Tissue
- Assay Buffer

1.3.2 Setup & Protocol

- First wash the flow chamber with 100 μ l of buffer, speed this up by holding a fiber-free tissue at the other side to 'pull' the buffer through. This movement is because of the capillary effect
- Add 30 μ l of BSA and let it incubate for 15 minutes
- Wash away residual BSA with 100 μ l assay buffer
- Add 30 μ l of your protein mix
- Close sides with nail polish to prevent water from evaporating