

Lab Notebook - Week 6 (17/7/17 - 23/7/17)

Project: NU iGEM 2017 Shared Project

Authors: Katerina Cheronis

Dates: 2017-07-17 to 2017-07-23

MONDAY, 7/17

Purpose: OMV purification from JC8031 (no antibiotic) and DeLisa ClyA-GFP (Cm+) fusion

OMVs were purified as described in the SBI protocol - 2 columns (one for GFP OMVs and ones for JC8031 no GFP)

1. Prepare clarified supernatant from bacterial culture
 - a. Culture the bacteria in its growth medium overnight at 37 deg.
 - b. Spin down bacteria at 5000xg for 15 minutes at 4 deg C
 - I. 4000xg for 17 minutes - centrifuge would not go any higher
 - c. Transfer the supernatant to a new flask and filter through a 0.45um filter
 - d. Spin down the supernatant at 5000xg for 15 minutes at 4 deg C
 - e. Transfer the supernatant to a new flask and filter through a 0.22um filter. The filtered supernatant is now ready for OMV isolation
2. Pack the column/bind OMVs
 - a. Pipette 200 mL of the resin onto the column
 - b. Equilibrate by adding 1mL of the Binding Buffer and allow the solution to flow through. Discard the flow through
 - c. Place the yellow cap onto the bottom of the column
 - d. Add 10mL of the clarified bacterial supernatant (prepared in step 1) to the resin and incubate on a rotating rack at 4 deg C for 3-4 hours to allow for OMV binding
 - I. Used centrifuge at a very low speed.
3. OMV Elution
 - a. Place the column onto a rack and allow the resin/supernatant to flow through (collect the flow through for analysis if desired)
 - b. Wash the resin with 10mL Binding Buffer 2 times. Discard the flow through.
 - c. Add 500uL Elution Buffer and collect in in 1.5ml Eppendorf tube
 - I. Repeat the elution step for a total of 5 times in separate tubes
4. Analysis
 - a. Perform downstream analysis of the five separate elutions (or pool if desired)

GEMSTONE Dye Protocol

Total time for completion: ~ 2 ++ hours

Useful for flow cytometry

- Resuspend 200 - 500ug protein equivalent of OMVs in 500uL PBS
- Add 1uL of the 500X or 2.5uL of the 250X labeling dye to the OMV preparation and incubate at 37degC with shaking for 20 minutes
- Add 167uL ExoQuick-TC to the solution and incubate at 4deg C for 2h-overnight
- Spin the eppendorf tube at 10000rpm for 10 minutes
- Carefully aspirate the supernatant from the corner of the tube
 - o No pellet formed.
- Resuspend the OMV pellet in 500 uL PBS and proceed with downstream applications - not done

Imaging: magnification not high enough, did not observe fluorescence or OMV aggregates

TUESDAY, 7/18

Preparation for Western blotting attempt #2

Day 1: 7/18/17

Retransform Competent JC8031:

Pos. Control for retransformation (10pg/ μ L DNA concentration): mRFP from iGEM distribution - Should appear **red** on cam plates

Neg. Control for retransformation: H₂O transformation - Should **not appear** on cam plates

Exp. Run for retransformation (10pg/μL DNA concentration): Purified plasmid DNA from Cas9-Dsba Miniprep

gRNA Design

Streak glycerol stocks of gRNA-mRFP and His6-Cas9

General:

- Made LB with Cm antibiotic
- Created Cm plates (approximately 40 plates)

Plate results:



WEDNESDAY, 7/19

General:

- Made LB : tetR & camR

First Retrasnformation Failed

Retry Retrasnformation

Experimental Setup:

Positive Control: Dh5a Competent Cells from NEB with Cas9-Dsba concentration 100 pg/uL

Experimental 1: JC8031 with with Cas9-Dsba concentration 100 pg/uL

Experimental 1: JC8031 with with Cas9-Dsba concentration 100 ng/uL

Experimental 1: JC8031 with with Cas9-Dsba concentration 100 pg/uL

ClyA-GFP Construct

Determined primers for sequencing ClyA-GFP - JBLF16 and CJK59

gRNA Design

Plates did not grow well, so no colonies were picked

Instead, directly made an 5mL overnight culture from glycerol stock of gRNA-mRFP and His6-Cas9

Prepared overnight cultures of:

mRFP - gRNA

His-Cas9 (Lulu)

His-Cas9 (Jack):

ClyA-GFP

DsbA - Cas9

OMV Experiment:

Retry Gemstone dyes to observe fluorescence:

Red dye used (Ruby)

OMVs taken from Tube 1: 2.953 mg/mL (Nanodrop 2000)

260/280 = 1.40

Tube 1 had the best protein concentration. 2nd 3rd and fourth elution concentrations were not as high.

Tube label: 29

Incubation time: approximately 4 hours

Stored in +4 for analysis

- In some papers OMVs are stored in -80
 - Where should they be stored?

THURSDAY, 7/20

Miniprep cultures prepared the previous day.

Miniprep Nanodrop results:

Using old Promega protocol

Attempt 1				
	A	B	C	D
1		Conc (ng/uL)	260/280	260/230
2	mRFP - gRNA	5.0	1.10	1.06
3	His-Cas9 Lulu	6.9	1.63	1.06
4	His-Cas9 Jack	40.0	1.64	0.71
5	ClyA-GFP	11.7	1.71	0.76
6	DsbA-Cas9	24.4	1.86	1.38

Using new protocol, Chelsea's revisions (new protocol available in miniprep kit):

Attempt 2				
	A	B	C	D
1		Conc (ng/uL)	260/280	260/230
2	mRFP - gRNA	9.4	1.71	0.66
3	His-Cas9 Lulu	44.8	1.91	2.21
4	His-Cas9 Jack	107.3	1.93	2.26
5	ClyA-GFP	27.8	1.94	2.60
6	DsbA-Cas9	168.4	1.90	2.22

Second miniprep attempt of ClyA - GFP was sent for sequencing (after 2:00 pm). Should receive early next week.

- Primers used:
 - CJGk59 (REV)// 5uL
 - JBLF16 (FW)//5uL
- Plasmid volume: 20uL (no dilution following nanodrop)

Prepared rich medium for Western blot with Jazzy (1/2 Dose)

Recipe for 1L 2xYT:

- Measure ~900ml of distilled H₂O.
- Add 16g Bacto Tryptone.
- Add 10g Bacto Yeast Extract.
- Add 5g NaCl.
- Adjust pH to 7.0 with 5N NaOH.
- Adjust to 1L with distilled H₂O.
- Sterilize by autoclaving.

Created new Tetracyclin 1000x stocks (10 stocks of 1 mL each) - based on protocol book concentrations

- In antibiotics cryo box (-20 deg). Aluminum foil because light sensitive.

Nanosight information from Devin:

Laser goes through the particles and uses Stokes Einstein equation for quantification

Size and concentration information

Sensitive to concentration - need to find best ratio (1×10^8 particles/mL)

- But we don't know OMV concentration following purification so this will require trial and error
- Western blot with flag tag would be the best for GFP detection
- Nanosight is able to detect fluorescence at certain excitations - worth a try
- Probably will try out using Nanosight early next week

gRNA Design:

Streatek plates of gRNA-mRFP from glycerol stock - Lulu

Prepared an overnight culture of mRFP-gRNA for miniprep the following day - Lulu

FRIDAY, 7/21

Plate results - Many aggregated colonies but there are some isolated ones to pick



Results from Retransformation 7/19/2017

Overall: Since colonies were left for more than 18 hours, the results are not to be used/trusted, will retry transformation again today. However, we will use the results as a reference to see what conditions to use for our retransformation

Positive Control: Dh5a Competent Cells from NEB with Cas9-Dsba concentration 100 pg/uL - **Contaminated with Ice, thrown out**

Experiment 1: JC8031 with with Cas9-Dsba concentration 100 pg/uL - **Few colonies**

Experiment 1: JC8031 with with Cas9-Dsba concentration 100 ng/uL - **Many colonies, grew extremely well**

Experiment 1: JC8031 with with Cas9-Dsba concentration 100 pg/uL - **Few colonies**

gRNA Experiment:

This experiment is to check the viability of the constructs **BBa_K2019002** and **BBa_K2019000**.

This follows the Biobricks submissions:

http://parts.igem.org/Part:BBa_K2019002

http://parts.igem.org/Part:BBa_K2019000

For mRFP-gRNA, the backbone is pSB1T3: http://beta.labgeni.us/registries/parts_registry/?part=BBa_K2019002

For His-Cas9, the backbone is pSB1C3: http://beta.labgeni.us/registries/parts_registry/?part=BBa_K2019000

Miniprep mRFP-gRNA cells from overnight culture (started on Thursday 7/20)

mRFP - gRNA				
	A	B	C	D
1	Amount of DNA during miniprep	DNA concentration	260/280	260/230
2	*High Concentration (#35)	349.1 ng/uL	1.89	2.08
3	**Low concentration (#36)	138.4 ng/uL	1.92	2.17

* Bigger pellet generated because of higher volume of bacteria culture used (not enough to use 2x1.5 mL for both preparations)

** 1.5 mL the first time and ~800 uL the second time the culture was centrifuged

Sequencing:

- Cas9-His was sent for sequencing with Forward (VF) and Reverse (VR) primer (10 uL template and 2 uL of each primer in separate tubes)

- gRNA - mRFP with the reverse primer only (VR)

Results for ClyA-GFP:

ClyA identified in sequence but GFP not there. Potentially downstream.

Next steps: sequence again with different primer.

Transformations of Cas9-DsbA construct from 2016 team

- 2 vials of Top10 comp. cells 2 vials of JC8031 comp cells
 - Thaw in ice for 15 mins
- 50 uL of each vial of comp cells was put in separate 1.7 mL microcentrifuge tubes
- 1 uL of plasmid 13 and plasmid 34 added to each cell strain (see bacterial database for more info)
- Vials were left in ice for 30 minutes
- Heat shock for 53 seconds followed by 5 minute incubation in ice
- Add 200 uL of SOC
- Aspirate and transfer to falcon tubes
- Incubate at 37 deg for 1 h and 15 m



SUNDAY, 7/23



Transformation results:

Note: the plates were left in the incubator for an additional day and then placed in the +4.

Started overnight cultures of JC8031 (no antibiotic) and JC8031 (ClyA-GFP) for induction with arabinose the following day (2 of each)