

## Lyophilization of Cas13a on paper

### Aim of the Experiment

In this experiment, the reaction mixture for Cas13a-dependent detection of RNA (described in protocol "plate reader experiments") is freeze-dried on paper to present a storable version of this assay.

### Materials

#### Preparation of Paper:

- Freeze-dryer
- Heat plate (magnetic stirrer, RH basic, IKA, Germany)
- Fluorescence-detector "Lightbringer" (designed by our team)
- Nuclease-free glass fiber filter paper (glass microfiber paper 934-AH RTU, Whatman, GE healthcare, Germany)
- Nuclease-free bovine serum albumine (BSA, 5%) (VWR life sciences, Germany)
- nuclease-free H<sub>2</sub>O (nf H<sub>2</sub>O, Sigma Aldrich, USA)

#### Cas13a reaction mixture without target RNA:

- Cas13a Lbu protein (from the species *Leptotrichia Buccalis*)
  - Lbu-specific and target-specific crRNA (DNA template for *in vitro* transcription: Biomers)
  - Target-RNA (extracted from cell lysate or *in vitro* transcribed)
  - Murine RNase inhibitor (NEB, Germany)
  - RNase A (IDT, Germany)
  - RNase processing buffer (10x, 200 mM HEPES [pH 6.8], 500 mM KCl, 50 mM MgCl<sub>2</sub>, 50 % glycerol)
  - nuclease-free H<sub>2</sub>O (nf H<sub>2</sub>O, Sigma Aldrich, USA)
-

Table 1: Cas13a mix final concentrations

Concentration	Chemicals
300 nM	Cas13a
50 nM	crRNA
1,5 U	RNase inhibitor
1x	RNase processing buffer

## Procedure

- Block filter paper (1 cm<sup>2</sup> per reaction) with BSA over night followed by thorough rinsing with nuclease-free water and dry on a heat plate at 80 °C for 20 min.
  - Add 15 µl Cas13a reaction mix without target RNA onto the blocked paper.
  - Immediately flash-freeze at -80 °C (alternatively: liquid nitrogen).
  - Place the loaded paper in the freeze-dryer and lyophilize at -60 °C for at least 4 h.
  - To perform an assay with this freeze-dried mixture, variable amounts of target-RNA-containing solution can be added and measured by the fluorescence detector.
-