

## PCR analysis of genomic DNA

Genomic DNA from *E. coli* MG1655 (Coli Genetic Stock Centre), DH5a (NEB) and Top10 (Invitrogen) was extracted using a Genomic DNA Purification Kit (ThermoFisher).

To confirm the presence/absence of the fimbriae operon, primers were designed to amplify fragments from the *fimA*, *fimD* and *fimH* coding sequences, Table 1.

CDS	Forward primer 5'-3'	Reverse primer 5'-3'	Product size
<i>fimA</i>	CTCTGGCAATCGTTGTTCTG	GAAGGTCGCATCCGCATTAG	518 bp
<i>fimD</i>	TCCACACTTCCCGATGAC	CACCGCTGACTCCGTAATAG	763 bp
<i>fimH</i>	GATGGGCTGGTCGGTAAATG	AGTCCCTACTGCTCCTAACG	763 bp

Table 1. List of primer sequences and the expected amplified product length.

Polymerase chain reactions (PCR) were set up as followed, using GoTaq DNA polymerase (Promega).

	1 rxn	Master mix for 4 rxn
5 x Green GoTaq reaction buffer	4 µl	12 µl
2 mM nucleotide mix	2 µl	8 µl
10 µM forward primer	2 µl	8 µl
10 µM reverse primer	2 µl	8 µl
5U/µl GoTaq DNA polymerase	0.1 µl	0.4 µl
10 ng/µl Genomic DNA	1 µl	
Nuclease free water	8.9 µl	35.6 µl

The PCR thermocycler was set up as follows:

Initial denaturation	95 °C	2 min
Denaturation	95 °C	30 s
Anneal	58 °C	30 s
Extension	72 °C	1 min
No. cycles	25	
Final extension	72 °C	5 min
Hold	10 °C	

PCR products were separated on a 1 % agarose gel and stained with SYBR Green for visualization.