

Time : 8 月 3 日

**1、Experimental purpose:** pGLO-Cas9 transformant colony PCR**2、Material**

Green Taq Mix (Vazyme)

Colony PCR primer: TCas-F: 5'- TCACACTTTGCTATGCCATAGC-3'

TCas-R: 5'- acggcatttagatagcgcacat-3'

**3、Experimental procedure**

(1) Four transformants were picked from the transformant plate for colony PCR verification;

(2) Dissolving the colonies in 20 ul of ddH<sub>2</sub>O, and picking a single colony containing the pGLO plasmid as a control, a total of 5 samples;

(3) Boiling the bacterial suspension in boiling water for 5 min, and the obtained solution is used as a template;

(4) Prepare the PCR system:


	Experimental group A1-A5/ ul	Control group B1-B5/ ul
2x Green Taq Mix	25	0
template	1	1
Primer F+R	2+2	2+2
ddH <sub>2</sub> O	20	45

1-4 is pGLO-Cas9 transformant and 5 is pGLO control

(5) PCR Program

94℃	5min	X30circle
94℃	30s	
50℃	30s	
72℃	4min20s	
72℃	5min	

(6) PCR product analysis by agarose gel electrophoresis

**4、Results**

The target product of the colony PCR was 3866 bp. The experimental group A1-A4 has about 4000 target strips, and the A5 no band indicates that the selected colonies are not false positive transformants containing pGLO. B1-

Time :   8  月   3  日

---

B4 has no strips, which excludes the effects of empty carrier running. The experimental results indicated that the selected single colonies were positive transformants containing pGLO-Cas9, and the next experiment will verify whether the protein expressed by the Cas9 gene is functional.