

WEDNESDAY, 11/07/2018

PCR Amplification of XylR gene from Top10 gDNA

Table1								
	A	B	C	D	E	F	G	H
1		Volume added (ul)						
2	Top10 gDNA	2.0			Initial denaturation	98°C	30s	
3	10 uM XylR_FP	2.5		34X	Denaturation	98°C	10s	
4	10 uM XylR_RP	2.5			Annealing	68°C	30s	
5	10 mM dNTP	1.0			Elongation	72°C	35s	1179bp
6	Q5 reaction buffer	5.0			Final elongation	72°C	120s	
7	Q5 polymerase	0.5				12°C	∞	
8	Nuclease-free water	36.5						
9	Total	50.0						

PCR Purification

Concentration of XylR : 90.3 ng/ul in 40 ul

Gel Electrophoresis

<image>

1+ kb band observed.

Restriction Digestion of pQE80L and XylR

Table2			
	A	B	C
1	Reaction	pQE80L (289.6 ng/ul)	XylR (90.3 ng/ul)
2		Volume added (ul)	
3	DNA	3.5 (~ 1 ug)	2.2 (~ 0.2 ug)
4	10X FastDigest Green Buffer	2.0	3.0
5	FastDigest BamHI	1.0	1.0
6	FastDigest HindIII	1.0	1.0
7	Nuclease-free water	12.5	22.8
8	Total	20	30

1. Add and mix the above.
2. Incubate at 37°C for 50 min.

**PCR Purification**

Concentration of pQE80L : 10.1 ng/ul in 40 ul

Concentration of XylR : 2.5 ng/ul in 40 ul

**Synchronized degradation tag characterization**

1. Transfer 200 µl of overnight culture into fresh 5 mL of LB+K (done the previous day).
2. Incubate at 37°C, 225 rpm for 1 hour (done the previous day).
3. Measure cell density at OD<sub>600</sub>.
4. Transfer 1ml of cell culture final OD<sub>600</sub> = 0.1, diluted using LB+K, into a 12-well plate. LB+K is added as the blank.
5. Repeat step 6 for triplicates.
6. Measure OD<sub>600</sub> and fluorescence intensity (535/600 nm) in a microplate reader using the following settings:  
 Temperature: Set-point 37°C  
 Plate In  
 Shake: Orbital for 10s  
 Read: (A) 600  
 Read: (F) 535,600  
 Excitation: 535  
 Emission: 600  
 Gain: 75  
 Read height: 6 mm  
 Plate Out
7. Cover the 12-well plate with a dark cloth.
8. Incubate at 25°C without shaking.
9. Repeat steps 8-10 every 1 hour over a period of 8 hours.

Experiment failed due to wrong incubation temperature and rpm used.

THURSDAY, 12/07/2018

**Ligation**

Table3					
	A	B	C	D	E
1	Reaction	Reaction 1: using previously cut pQE80L (52.4 ng/ul)	Reaction 2: using cut pQE80L from above (10.1 ng/ul)		
2		Volume added (ul)			
3	pQE80L	1.0	5.0	~ 37.5 ng	3:1 ratio, using NEBioCalculator for ligation
4	XylR	15.0	15.0	~ 50 ng	
5	T4 DNA ligase	1.0	1.0		
6	T4 DNA ligase buffer	2.0	2.5		
7	Nuclease-free water	1.0	1.5		
8	Total	20.0	25.0		

Incubation at 25°C for 35 min.

**Transformation into DH5α**

Transform 10 ul from reaction 1 and 15 ul from reaction 2 separately to competent DH5α.

Spread onto AmpR plates: DH5α-pQE80L-XylR-1 and DH5α-pQE80L-XylR-2.

FRIDAY, 13/07/2018

**Colony patching**

Satellite colonies observed.

Patch 10 colonies from both DH5α-pQE80L-XylR-1 and DH5α-pQE80L-XylR-2 plates onto new agar plates.

Incubate at 37°C.

**Colony PCR**

Table4								
	A	B	C	D	E	F	G	H
1		Volume added (ul)						
2	10 uM XylR_FP	1.0			Initial denaturation	95°C	10s	
3	10 uM XylR_RP	1.0		34X	Denaturation	95°C	10s	
4	10 mM dNTP	1.0			Annealing	<b>55°C</b>	30s	
5	10X Std Taq Buffer	5.0			Elongation	72°C	<b>72s</b>	1179bp
6	Taq polymerase	0.5			Final elongation	72°C	300s	
7	Nuclease-free water	41.5				12°C	∞	
8	Total	50.0						

1. Prepare two of the above reaction.
2. Mix and transfer 10 ul into four other tubes.
3. Pick colonies 1 to 5 from DH5α-pQE80L-XylR-1 and DH5α-pQE80L-XylR-2 respectively and transfer to respective tubes.

**Gel Electrophoresis**

1kb band observed for all lanes except lane 10 (DH5α-pQE80L-XylR-2 colony 5).

Intensity of band highest for lane 3 (DH5α-pQE80L-XylR-1 colony 3), followed by lane 5 (DH5α-pQE80L-XylR-1 colony 5).

**Quadrant Streak** DH5α-pQE80L-XylR-1 colony 3 onto new AmpR agar plate and incubate at 37°C overnight.

**Synchronized degradation tag characterization**

1. Transfer 200 µl of overnight culture into fresh 5 mL of LB+K (done the previous day).
2. Incubate at 37°C, 225 rpm for 1 hour (done the previous day).
3. Measure cell density at OD<sub>600</sub>.
4. Transfer 1ml of cell culture final OD<sub>600</sub> = 0.1, diluted using LB+K, into a 12-well plate. LB+K is added as the blank.
5. Repeat step 6 for triplicates.
6. Measure OD<sub>600</sub> and fluorescence intensity (535/600 nm) in a microplate reader using the following settings:  
 Temperature: Set-point 37°C  
 Plate In  
 Shake: Orbital for 10s  
 Read: (A) 600  
 Read: (F) 535,600  
 Excitation: 535  
 Emission: 600  
 Gain: 75

Read height: 6 mm

Plate Out

7. Cover the 12-well plate with a dark cloth.
8. Incubate at 37°C, 120 rpm.
9. Repeat steps 8-10 every 1 hour over a period of 8 hours.

---

MONDAY, 16/07/2018

**Inoculation of DH5α-pQE80L-XylR** into 10 ml of LB broth and 10 ul of ampicillin

---

TUESDAY, 17/07/2018

**Preparation of DH5α-pQE80L-XylR glycerol stock**

**Plasmid Extraction**

Concentration of pQE80L-XylR : 132.7 ng/ul in 40 ul

To be sent for sequencing.

**Synchronized degradation tag characterization (re-run)**

1. Transfer 200 µl of overnight culture into fresh 5 mL of LB+K (done the previous day).
2. Incubate at 37°C, 225 rpm for 1 hour (done the previous day).
3. Measure cell density at OD<sub>600</sub>.
4. Transfer 1ml of cell culture final OD<sub>600</sub> = 0.1, diluted using LB+K, into a 12-well plate. LB+K is added as the blank.
5. Repeat step 6 for triplicates.
6. Measure OD<sub>600</sub> and fluorescence intensity (535/600 nm) in a microplate reader using the following settings:

Temperature: Set-point 37°C

Plate In

Shake: Orbital for 10s

Read: (A) 600

Read: (F) 535,600

Excitation: 535

Emission: 600

Gain: 75

Read height: 6 mm

Plate Out

7. Cover the 12-well plate with a dark cloth.
8. Incubate at 25°C without shaking.
9. Repeat steps 8-10 every 1 hour over a period of 8 hours.

---

THURSDAY, 19/07/2018

XylR sequence confirmed.

**PCR Mutagenesis**

Table5

	A	B	C	D	E	F	G	H	I	J
1		Reaction 1 (R121C)	Reaction 2 (P363S)	Reaction 3 (R121C, P363S)						
2		Volume added (ul)								
3	pQE80L-XyIR (132.7 ng/ul)	0.5	0.5	0.5			Initial denaturation	95°C	120s	
4	10 uM R121C_FP	1.0	0	1.0		30X	Denaturation	95°C	20s	
5	10 uM R121C_RP	1.0	0	1.0			Annealing	<b>58°C</b>	20s	
6	10 uM P363S_FP	0	1.0	1.0			Elongation	72°C	<b>180s</b>	~6kb
7	10 uM P363S_RP	0	1.0	1.0			Final elongation	72°C	300s	
8	2.5 mM dNTP	4.0	4.0	4.0				12°C	∞	
9	5X TranStart Fast Pfu Buffer	10.0	10.0	10.0						
10	TranStart Fast Pfu DNA polymerase	1.0	1.0	1.0						
11	Nuclease-free water	32.5	32.5	30.5						
12	Total	50.0	50.0	50.0						

**DpnI Digestion**

Table7

	A	B
1		Volume added (ul)
2	DNA from PCR mutagenesis	30.0
3	10X CutSmart Buffer	3.5
4	DpnI	1.0
5	Nuclease-free water	0.5
6	Total	35.0

The above was conducted for all three reactions. Incubate at 37°C for 1 hr.

TUESDAY, 24/07/2018

**PCR Purification**

Concentration of DH5α-pQE80L-XyIR\*-1 (R121C) : 112.6 ng/ul in 40 ul

Concentration of DH5α-pQE80L-XyIR\*-2 (P363S) : 57.1 ng/ul in 40 ul

Concentration of DH5α-pQE80L-XyIR\*-3 (R121C, P363S) : 100.8 ng/ul in 40 ul

**Transformation into DH5α**

2 ul of DNA added to competent DH5α.

DH5α-pQE80L-XyIR\*-1 (R121C)

DH5α-pQE80L-XyIR\*-2 (P363S)

DH5α-pQE80L-XyIR\*-3 (R121C, P363S)

WEDNESDAY, 25/07/2018

Colony Patching and Inoculation

DH5α-pQE80L-XylR\*-1 (R121C) : many colonies  
DH5α-pQE80L-XylR\*-2 (P363S) : no colony  
DH5α-pQE80L-XylR\*-3 (R121C, P363S) : 7 colonies

Patch 7 colonies from DH5α-pQE80L-XylR\*-1 and -3 onto new plate.  
Incubate at 37°C.  
Inoculate colony 1 from both DH5α-pQE80L-XylR\*-1 and -3 in 10 ml of LB broth with 10 ul of ampicillin.

THURSDAY, 26/07/2018

Plasmid Extraction

DH5α-pQE80L-XylR\*-1 colony 1 : 118.1 ng/ul in 40 ul  
DH5α-pQE80L-XylR\*-3 colony 1 : 85.3 ng/ul in 40 ul  
To be sent for sequencing

Inoculation

Inoculate colonies 6 and 7 from both DH5α-pQE80L-XylR\*-1, and colonies 4 and 6 from DH5α-pQE80L-XylR\*-3 in 10 ml of LB broth with 10 ul of ampicillin.

FRIDAY, 27/07/2018

Plasmid Extraction

DH5α-pQE80L-XylR\*-1 colony 6 : 102.8 ng/ul in 40ul  
DH5α-pQE80L-XylR\*-1 colony 7 : 97.8 ng/ul in 40 ul  
DH5α-pQE80L-XylR\*-3 colony 4 : 119.4 ng/ul in 40 ul  
DH5α-pQE80L-XylR\*-3 colony 6 : 77.4 ng/ul in 40 ul  
To be sent for sequencing.

TUESDAY, 31/07/2018

P363S Mutagenesis

Table6									
	A	B	C	D	E	F	G	H	I
1		Reaction 1 (R121C)	Reaction 2 (P363S)						
2		Volume added (ul)							
3	pQE80L-XylR*-1 colony 1 (118.1 ng/ul)	0.5	0.5						
4	10 uM P363S_FP	1.0	10.0			Initial denaturation	95°C	120s	
5	10 uM P363S_RP	1.0	10.0		30X	Denaturation	95°C	20s	
6	2.5 mM dNTP	4.0	4.0			Annealing	55°C	20s	
7	5X TranStart Fast Pfu Buffer	10.0	10.0			Elongation	72°C	180s	~6kb
8	TranStart Fast Pfu DNA polymerase	1.0	1.0			Final elongation	72°C	300s	
9	Nuclease-free water	32.5	14.5				12°C	∞	
10	Total	50.0	50.0						

Mutagenesis for P363S re-attempted using higher primer concentrations (Reaction 2) and lower Tm (55°C).

