

PCR protocol

1. assemble all reaction components.

component	50µl reaction
Buffer(2×)	25µl
DNA polymerase	1µl
Forward primer(20µM)	1µl
Reverse primer(20µM)	1µl
dNTP	1µl
Template DNA	1µl
ddH ₂ O	20µl

Note:

(1) stick pipette into liquid while adding components.

(2) get the liquid on the tube wall down by centrifuge to decrease the loss of fluid.

(3) have total volume 1~2µl lesser than 50µl in final components assembly because you will lose some fluid while transferring.

2. PCR reaction.

step	temperature	time
Initial denaturation	95°C	5min
Cycle reaction(30 times)	95°C	30s
	56°C	30s
	72°C	1kb/min
Final extension	72°C	5min
hold	4°C	/