

Genome editing via the CRISPR-Cas9 system

Genome Editing

- Generating *E. coli* BL21 DE3 harboring pCas
 - Do a heatshock transformation of pCas with *E. coli* BL21 DE3
 - Grow the culture overnight in LB media and kanamycin at 30 °C.
 - Generate electrocompetent *E. coli* BL21 DE3 cells harboring pCas.
- Add Arabinose (10 mM) to the culture for lambda-Red induction.
- Transformation via electroporation of pTarget
 - Thaw 50 µL electrocompetent *E. coli* cells on ice.
 - Add 100 ng pTarget and 400 ng of the donor DNA (PCR fragment).
 - Store cells on ice for 1 minute.
 - Electroporate at $U = 2.5 \text{ kV}$, $C = 25 \text{ µF}$, $R = 400 \text{ Ω}$
 - Recover cells at 30 °C for 1 hour.
 - Plate on selective LB agar containing kanamycin (50 mg L⁻¹) and spectinomycin (50 mg L⁻¹).
 - Incubate overnight at 30 °C.
 - Verify positive transformants by colony PCR and DNA sequencing.

Plasmid curing

- Inoculate a positive colony, harboring pCas and pTarget in 2 mL LB media containing kanamycin (50 mg L⁻¹) and IPTG (0.5 mM).
- Incubate at 30 °C for 8 to 16 hours.
- Plate on selective LB agar containing kanamycin (50 mg L⁻¹)
- To cure pCas, grow the cells overnight at 37 °C non-selectively.

From: iGEM Bielefeld-CeBiTec