



#### **QIAGEN- Miniprep:**

- 1) Pellet 5 ml bacterial overnight culture by centrifugation at 4700 x g for 3 minutes at room temperature.
- 2) Resuspend pelleted bacterial cells in 250 µl Buffer P1 (Resuspension Buffer) and transfer to a microcentrifuge tube.
- 3) Add 250 µl Buffer P2 (Lysis Buffer) and mix thoroughly by inverting the tube until the solution becomes clear.
- 4) Add 350 µl of Buffer N3 (Neutralization Buffer) and mix immediately and thoroughly by inverting the tube 4-6 times.
- 5) Centrifuge for 10 minutes at 13000 rpm in a table-top microcentrifuge.
- 6) Apply 800 µl supernatant from step 5 to the QIAprep 2.0 spin column by pipetting. Centrifuge for 30-60 seconds and discard the flow-through.
- 7) Wash the column by adding 0.5 ml Buffer PB (Binding Buffer). Centrifuge 30-60 seconds and discard the flow-through.
- 8) Wash the column by adding 0.75 ml Buffer PE (Washing Buffer). Centrifuge for 30-60 second and discard the flow-through. Transfer the column to the collection tube.
- 9) Centrifuge for 1 min to remove residual wash buffer.
- 10) Place column in a clean 1.5 ml microcentrifuge tube. Add 50 µl ddH<sub>2</sub>O to the center of the column, let stand for 1 min, and centrifuge for 1 min.
- 11) Measure concentration of the sample with Nanodrop.
- 12) Store sample at -20°C.