

Based on the results we have achieved, we chose lacI-P_{lac} to carry out the next step. In this system, the arabinose-controlled repressor protein is utilized to regulate the transcription of sgRNA.

2018-8-21

One step cloning and transformation

Materials (20 µl Reaction) :

5 × CE Multis Buffer	4µl
pGLO-lacI	2µl
sgRNA	1µl
Exnase Multis	2µl
ddH ₂ O	11µl

2018-8-23

Colony PCR

Materials (20 µl Reaction) :

Green Taq Mix	25µl
P1	2µl
P2	2µl
Boiled bacteria solution	1µl
ddH ₂ O	20µl

PCR condition		pGLO-lacI-sgRNA
Step one		94°C 300s
30 cycles	Step two	94°C 30s
	Step three	48°C 30s
	Step four	72°C 60s
Step five		72°C 420s

Primers:

P1 5'-ATACGACGATACCGAAGACA-3'

P2 5'-CAAGTTGATAACGGACTAGC-3'

Here is the result and the correct is to get a 633bp fragment.

