


# Arabinose Promoter (BBa\_I0500) + Agr C & A

## Plasmid Generation Protocol



1. Miniprep overnight cultures of **BBa\_I0500** (Arabinose Promoter) and **BBa\_I746107** (Agr C & A). This is your backbone and insert.
2. Digest BBa\_I0500 plasmid with **SpeI** and **PstI** at **37°C** for **60mins**. Use the following mixture:

Molecular Grade Water	Calculated. Use to make reaction total volume = 50 $\mu$ L
10x NEB Buffer 4	5 $\mu$ L
DNA	1 $\mu$ g
SpeI	0.5 $\mu$ L
PstI	0.5 $\mu$ L

3. Digest BBa\_I746107 (Agr C & A) plasmid with **XbaI** and **PstI** at **37°C** for **60mins**. Use the following mixture:

Molecular Grade Water	Calculated. Use to make reaction total volume = 50 $\mu$ L
10x NEB Buffer 4	5 $\mu$ L
DNA	1 $\mu$ g
XbaI	0.5 $\mu$ L
PstI	0.5 $\mu$ L

4. Clean up all DNA fragments with PCR clean-up kit.
5. Using the “Ligation Template” excel sheet, calculate the amount of each component to combine in a ligation mix for an Insert:Backbone ratio of **4:1** with total DNA concentration of **7.5-10.0ng/uL**.

Incubate at 25°C for 1 hr. Heat inactivate at 80 C for 10 mins. Place on ice.

6. **Transformation** (Need 2x LB+Cm20 plates. Check for plates. Check for hockey spreaders/glass beads):

- Set water bath to 42°C.

- Remove LB+Cm20 plates from 4°C and allow them to come to RT.
  - Thaw chemically competent cells on ice. Leave in microcentrifuge tube.
  - Add **5µL** of each ligation mix to individual chemically competent cells.
  - Incubate on ice for 30 min.
  - Heat shock cells for 60 sec at 42°C without shaking.
  - Place on ice for 2 minutes.
  - Aseptically (by the fire or in the hood) add **250µL** of SOC media to the tube. Cap tightly.
  - Place tube horizontally in shaker. Incubate at 37°C and 225 rpm for 1 to 2 hr.
  - In the laminar hood, spread 100uL of transformants onto LB+ Cm20 plates.
  - Leave plates in 37°C incubator overnight. Store remaining liquid cultures in 4°C.
7. Next day, pick and annotate (give each colony a #) **16** colonies onto new LB+ Cm20 plates. Place new plates in 37°C overnight.
8. **Verify all 16 colonies** via PCR with **Taq DNA polymerase** using these primers: VF2 and VR. Use following mixture:

Component	50µl reaction
Molecular grade H2O	Added first. 38.5µL
10x Standard Taq Buffer	5µl
10 mM dNTPs	1µl
Sequencing VF2 (10µM)	2.5µl
Sequencing VR (10µM)	2.5µl
<i>Taq DNA polymerase</i>	Added last. 0.5µl

Sequencing VF2: TGCCACCTGACGTCTAAGAA [Tm = 55.9]  
 Sequencing VR: ATTACCGCCTTTGAGTGAGC [Tm = 55.3]

Use the following thermalcycler protocol:

Cycle Step	3-step protocol		Cycles
	Temp (°C)	Time	
Initial denaturation	95	5 min	1
Denaturation	95	30s	<b>35</b>
Annealing	<b>68</b>	30s	
Extension	68	<i>4min 30sec</i>	
Final Extension	68	10min	1
Hold	4	hold	

9. Run **10µL** of all 16 PCR products onto gel (110V for 60mins). [Do you see the expected size(s)? expected size is 1109 bp]
10. **If** you find colonies with the aforementioned sizes, **PICK BEST TWO MUTANTS** and start overnight cultures in 2mL LB + Kan50. Leave them to grow in the 37°C incubator overnight.
11. Next morning; miniprep them. Measure their concentrations with the nanodrop. Send them for sequencing.