

Week 14

Vienna RNA thermoswitches

AIM:

Make pre-made RNA thermoswitches for GFP nano (igA) constructs.

REAGENTS USED:

Tab. 1 List of reagents used in the experiment

Name of reagent
100pb Plus DNA ladder (Thermo Scientific)
GeneJet PCR purification Kit
GeneJet Gel Purification Kit
DNA loading dye 6x
Isopropanol
DH5 α competent cells
LB media
LB + agar plates
Chloramphenicol
SOC

Tab. 2 List of primers used in the experiment

Name of primer pair	Forward sequence	Reverse sequence
GJ2	GAATTCGCGGCCGCTTCTAGAGGTAT CCCTCACTTACTAG	CTGCAGCGGCCGCTACTAGTAGCTAGCGGGTATA TCTCCTTCTG
GJ3	GAATTCGCGGCCGCTTCTAGAGGGTC CCTCACTTACT	CTGCAGCGGCCGCTACTAGTAGCTAGCGGGTATA TCTCCTTCTG
GJ6	GAATTCGCGGCCGCTTCTAGAGTCTC CTTCACTAGTAA	CTGCAGCGGCCGCTACTAGTAGCTAGCGGGTATA TCTCCTTTT
GJ9	GAATTCGCGGCCGCTTCTAGAGCTCC T TACTAGTCTGC	CTGCAGCGGCCGCTACTAGTAGCTAGCGGGTATA TCTCCTTCTG
GJ10	GAATTCGCGGCCGCTTCTAGAGTCTC CTTACTAGTCTGC	CTGCAGCGGCCGCTACTAGTAGCTAGCGGGTATA TCTCCTTCTG
Sw2Core	AAGAGAAAGCGGTCCGTCTAAGGAG CCTGAGGT	ACCTCAGGCTCCTTAGACGGACCGCTTTCTCTT
Sw2BB	GAATTCGCGGCCGCTTCTAGAGAAG AGAAAGCGGTCCG	CTGCAGCGGCCGCTACTAGTAACCTCAGGCTCCTT AG
Sw3Core	AGTGCGGTGTATTTATAACAAGGAGC GGGATGG	CCATCCCGCTCCTTGTATATAAATACACCGCACT
Sw3BB	GAATTCGCGGCCGCTTCTAGAGAGTG CGGTGTATTTAT	CTGCAGCGGCCGCTACTAGTACCATCCCGCTCCT
Sw6Core	GTGCGGTGCCGTGGTCACTAAGGAG CCCTTGGG	CCCAAGGGCTCCTTAGTGACCACGGCACC GCAC
Sw6BB	GAATTCGCGGCCGCTTCTAGAGGTGC GGTGCCGT	CTGCAGCGGCCGCTACTAGTACCCAAGGGCTCCT
Sw7Core	GATCGCCCCTAGGAGGAACAAGGAG CTATCGTT	AACGATAGCTCCTTGTTCCTCCTAGGGGCGATC

Sw7BB	GAATTCGCGGCCGCTTCTAGAGGATC GCCCCTAGGA	CTGCAGCGGCCGCTACTAGTAAACGATAGCTCCTT GT
Sw9Core	GAGGTTTTATAGTCCTTCTAAGGAGA TTGAGCT	AGCTCAATCTCCTTAGAAGGACTATAAACCTC
Sw9BB	GAATTCGCGGCCGCTTCTAGAGGAG GTTTTATAGTCCTTC	CTGCAGCGGCCGCTACTAGTAAGCTCAATCTCCTT AG
Sw11Core	GGTGCGGTGTATAAGGTGAAAGGAG TCATTGGG	CCCAATGACTCCTTTCACCTTATACACCGCACC
Sw11BB	GAATTCGCGGCCGCTTCTAGAGGGT GCGGTGTATAAG	CTGCAGCGGCCGCTACTAGTACCCAATGACTCCTT TC

EXPERIMENT DESCRIPTION:

Due the lack of time, RNA thermoswitches were synthesized using 4 primers. Core primers make the thermoswitch sequence, and BB primers insert biobrick sites at the beginning and the end of the gene.

EXPERIMENT PROTOCOL:

- Run PCR using for primers: **BBforward**, **BBreverse**, **Coreforward** and **Corereverse** primers to synthesize RNA thermoswitches and use GJx primers, to insert biobrick sites on the old RNA thermoswitches.
- PCR reaction

Reagent	Amount

Super Fi MM x2	25 μ L
Rcore+Fcore / GJx forward	0.5+0.5 μ L (2.5 nM each) / 2.5 (10 nM)
BBpre+BBsuff / GJx reverse	2.5+2.5 μ L (10 nM each) / 2.5 (10 nM)
H ₂ O	19 / 20 μ L
Total:	50 μ L

- Run gel electrophoresis.
- Purify PCR products using GeneJet PCR Purification Kit (ThermoFisher scientific).
- Digest GJ_x with XbaI and NheI, Sw_x with XbaI and BclI restriction enzymes and pSB1C3 plasmid with XbaI and BclI using FastAP. Digest for 1 hour.
- Digestion reaction

Reagent	Amount
DNA	200 ng / 3000 ng
FastAP	- / 5 μ L
10x Fast Digest buffer	2 / 5 μ L
XbaI	0.5 / 5 μ L
BclI (NheI for Sw _x)	0.5 / 5 μ L
H ₂ O	To 20 / 50 μ L
Total:	20 / 50 μ L

- Run gel electrophoresis and purify cut plasmid using the GeneJet Gel Purification Kit (ThermoFisher scientific).
- Ligate RNA thermoswitches into cut plasmid for 1 hour.
- Ligation reaction

Reagent	Amount
5x Rapid ligation buffer	4 μ L

T4 ligase	1 μ L
DNA	60 ng
Insert	Xng (7:1 compared to DNA by molar mas)
H ₂ O	Up to 20 μ L
Total:	20 μ L

10. Transform DH5 α competent cells with ligation mix.
11. Plate cells on LB+ agar plates with chloramphenicol.
12. Grow in 37°C for 16 h.
13. Run colony PCR using BBforward and BBreverse primers for Sw_x and for GJx use forward and reverse primers.

Reagent	Amount
Primer BBreverse	0.5 μ L
Primer BBforward	0.5 μ L
DreamTaq MM x2	2 μ L
H ₂ O	4 μ L
Total:	10 μ L

14. Two step PCR.

1 cycle	95°C	2 minutes
6 cycles	95°C	15 seconds
	50/54/57°C	15 seconds

	72°C	10 seconds
1 cycle	72°C	3 minutes
30 cycle	95°C	15 seconds
	72°C	10 seconds
1 cycle	72°C	5 minutes
1 cycle	4°C	indefinite period

15. Data was analysed by running the electrophoresis gel after PCR reaction, digestion and colony PCR.

RESULTS:

After PCR reaction products of 88 bp- in first line, 87 bp-second line, 102 bp-third line, 81 bp- forth line, 82 bp-fifth line and 76 bp in the rest of lines were expected. You can see them in Figure 1.




Fig. 1 Electrophoresis gel of PCR. 1- GJ2, 2- GJ3, 3- GJ6, 4- GJ9, 5- GJ10, 6- Sw2, 7- Sw3, 8- Sw6, 9- Sw7, 10- Sw9, 11- Sw11

After PCR, products were cut, purified and ligated to pSB1C3 plasmid. After transformation and 16 h growth in 37°C, transformed bacteria were checked by running colony PCR. All colonies came out to be positive, except a few ones with non-specific PCR products (those colonies were not transfer to mini preps for plasmid multiplication). You can see the result in a picture two. PCR product length is the same as in the first PCR.




Fig. 2 Colony PCR gel electrophoresis.

CONCLUSIONS:

1. Designed primers were specific for RNA thermostwitches production, as well as for the colony PCR.
2. Both ligation and the bacterial transformation showed the plasmid with genes of interest.