

# Strep-Tag protein purification

The Strep-Tag II is a protein tag that enables the purification and detection of recombinant proteins by affinity chromatography. It is a synthetic peptide which consists of 8 amino acids and can be expressed N- or C-terminal with a fusion protein. Strep-tag II has a strong affinity for strep-tactin, which is used for the purification of fusion proteins via affinity chromatography columns.<sup>1</sup>

## Materials

Buffer A Strep-Tag (pH: 7.4; 20mM NaP; 280mM NaCl; 6mM KCl)

Buffer B Strep-Tag (pH: 7.4; 20mM NaP; 280mM NaCl; 6mM KCl; 2.5mM Desthiobiotin)

Lysozyme

DNase

Benzamidine

PMSF

MilliQ Water

Dialysis Buffer (DIA)

## Instruments:

Sonifier

Centrifuge

Liquid pump

ÄKTApurifier

Glass flask

Measuring cylinder

Stir bar

Dialysis Equipment

## Procedure

Lysis

Add 3-5 times as much mL buffer A Strep-Tag as mg pellet.

Add a spatula tip lysozyme

Add DNase: 1/10000 of lysate volume

Add benzamidine: 5mM final concentration

Add PMSF: 0.5mM final concentration

Sonification:

Amplitude 50% (small peak); 30% (large peak)

5:00 min

0,5 Pulse on

2,0 Pulse off

Remove lysate -> take sample "Lysis Sample"

Centrifuge lysate 17k 50min 4°C

supernatant in glass flask -> take sample "supernatant sample"

resuspend pellet -> take sample "pellet sample"

Equilibrate Strep column (at 2 mL/min)

2-3x volume of column wash with MilliQ water

2-3x volume of column wash with buffer A Strep-Tag

circulate 45min lysate through column -> take sample "flow through"

ÄKTA

Connect buffer A and B Strep-tag at ÄKTA

Start PumpWash

Connecting the column to ÄKTA

Run buffer A through column (until UV value has reached basal line)

Add 100% buffer B (reduce fraction size)

Removing fractions

Start PumpWash

Nanodrop measurement of fractions

Determination of the concentration via the specific extinction coefficient

(<https://web.expasy.org/protparam/>)

Dialysis of samples

Add 1L dialysis buffer (DIA) to the measuring cylinder.

Moisten dialysisbag 4-7kDa with DIA buffer.

Attach knots and brackets and collected fractions as follows

Allow to stand overnight at 4 °C while stirring.




Figure 1: Dialysis of samples

New concentration determination with the Nanodrop

Freeze samples with liquid nitrogen and store at -80°C (add glycerol, if necessary)

## Trouble shooting

The column should be thoroughly and long washed in the ÄKTA.

The dialysis bag must be in a moving liquid to wash the buffer B out.

Fast working is important in this process to prevent additional denaturation of the proteins.

The work should take place from the lysis at 4°C.

## References

[1] A Skerra: Das Strep-tag als molekulares Werkzeug zur Hochdurchsatz-Proteinreinigung in der Proteomforschung. In: BIOSpektrum. 2, 2003, S. 189-192.

[2] Purification protocol provided by AG Binz

## Figures:

[1]: <http://de.spectrumlabs.com/dialysis/biotech.html>