

Protocol Name: Glycerol Stock Preserves

Category: General

Date: 19/09/2018

Author: Connor Trotter

Source(s): Addgene (2016) Creating Bacterial Glycerol Stocks For Long-Term Storage of Plasmids
<https://www.addgene.org/protocols/create-glycerol-stock/> (Accessed: 01/07/2018)

Time Required: 1 hour preparation day before, 30 minutes on the day

Additional Notes: Timing may be dependent upon growth rate of bacteria

Materials:

15ml LB broth (per species)

1x 50ml falcon tube (per species)

10ml 50% glycerol

Sterile 1.5ml Eppendorf tubes

Liquid Nitrogen (Optional)

100-1000µl pipette and appropriate sterile tips

Procedure:

- 1) The day prior to producing the stock preserves, inoculate 15ml LB in a 50ml falcon tube with a single colony of the desired species from an agar plate.
 - a. Incubate at appropriate temperature and rotations per minute for 12-16 hours
- 2) Observe falcon tube to determine if there has been substantial bacterial growth
- 3) Rock the falcon tube to resuspend any sediment cells
- 4) Add 500µl of overnight culture to 10 (or less) 1.5ml Eppendorf tubes
- 5) Add 500µl of 50% glycerol to the Eppendorf tubes containing the overnight culture
- 6) (Optional, dependent upon bacteria) Taking appropriate safety measures, freeze the Eppendorf tubes in liquid nitrogen for 10 seconds then retrieve
- 7) Store the Eppendorf tubes at -80°C
- 8) To retrieve bacteria from a stock, use a sterile loop or pipette tip to scrape some bacteria off the top. The object can then be used to plate bacteria or inoculate liquid culture as normal
 - a. Aim to not thaw the stock as this will reduce the lifespan of the stock if refrozen