

Time : 8 月 2 日

1、 Experimental purpose: pTF-p15A transformant colony PCR

2、 Material

Green Taq Mix (Vazyme)

Primer:

F: 5'— cacgatcgacattgatctgg—3'

R: 5'— GTGCATTCCGCTGTTATGG—3'

3、 Experimental procedure

(1) Pick 8 transformants for colony PCR verification;

(2) After picking the selected single colonies and spotting on another spectinomycin-resistant plate, the colonies were dissolved in 20 ul of ddH₂O, and a single colony containing the pTargetF plasmid was picked as a control, a total of 9 samples;

(3) The bacterial suspension was boiled in boiling water for 5 min, and the resulting solution was used as a template;

(4) Prepare the PCR system:


	Experimental group A1-A9/ ul	Control group B1-B9/ ul
2x Green Taq Mix	25	0
template	1	1
Perimer F+R	2+2	2+2
ddH ₂ O	20	45

(5)PCR program

94℃	5min	X30 circle
94℃	30s	
50℃	30s	
72℃	70s	
72℃	5min	

(6) PCR product analysis by agarose gel electrophoresis

4、 Results



Time : 8 月 2 日

The target product of colony PCR was 1033 bp. The electrophoresis results showed that the PCR products of the experimental group 1-8 were slightly larger than 1000 bp, and there was no band in the control group. There were no bands in the pTargetF control group A9 and B9, so 1-8 were non-empty carriers. This can be used to preliminarily demonstrate that pTF-p15A was successfully constructed. The next step will be to verify the function.