



New England BioLabd- Phusion PCR:

1) Set up the following reaction on ice:

COMPONENT	AMOUNT
5X Phusion HF Buffer	10 μ l
10mM dNTPs	1 μ l
10mM Forward Primer	2,5 μ l
10mM Reverse Primer	2,5 μ l
Phusion DNA Polymerase	0,5 μ l
Template DNA	0,5 μ l
ddH ₂ O	33 μ l
Total amount	50 μl

2) Split 50 μ l into two tubes with 25 μ l.

3) Set up following PCR program on thermocycler:

STEP	TEMP	TIME
Initial Denaturation	98°C	30 seconds
25-30 cycles	98°C 68°C 72°C	10 seconds 10 seconds 15-30 seconds/kb
Final Extension	72°C	10 seconds
Hold	12°C	