

The experiment design of 2018 OUC-China collaboration

1. PURPOSE

Measuring the fluorescence and OD600 of the target strain

2. Materials

- a. LB medium (peptone 10g/L, yeast extract 5g/L, NaCl 10g/L)
- b. M9 medium (12.8 g/L Na₂HPO₄·4.7H₂O, 3g/L KH₂PO₄, 0.5 g/L NaCl, 1.67g/L NH₄Cl, 1 mM thiamine hydrochloride, 0.4% glucose, 2 mM MgSO₄, 0.1 mM CaCl₂)
- c. IPTG (0.125mM)
- d. 96-well plates (the picture below)
- e. target strain and DH5 α



3. Protocol

- a. Take the glycerine bacteria plate (30-40 μ l is enough) and culture overnight.
- b. Single colony is selected and cultured overnight in 50ml LB medium (150rpm, 37°C).
- c. Take 1ml LB bacteria solution into 99ml M9 medium.
- d. Add IPTG after 4h culture (200rpm, 37°C).
- e. Starting with the addition of IPTG, OD600 and fluorescence were measured every hour, and 8 sets of data were measured.

4. Attentions:

- a. The bacteria added with IPTG, without IPTG, DH5, M9 medium should be measured.
- b. At the time of measurement, 200 μ l was taken from the bottle each time in the perforated plate, and 3 repeated groups were taken each time for each bottle bacteria.
- c. When selecting single colonies, if the bacteria solution is too thick and can be diluted appropriately, three single colonies should be selected for biological repetition.
- d. Fluorescence (485nm excitation light and 520nm emission light) and OD600.