

Protocol Name: Cell Counting

Category: Chemotaxis

Date: 01/09/2018

Author: Connor Trotter

Source(s):

Time Required: 1 day preparation + 1 day Experimentation

Additional Notes:

Materials

1x Glass Haemocytometer

15ml LB Media

1x 50ml Falcon Tube

Sterile Inoculation Loops

Sterile Pipette Tips

Coverslips

5% Chemgene/Disinfectant

70% Ethanol

Procedure

- 1) In sterile conditions, add 15ml LB media to a 50ml falcon tube
- 2) Using a sterile inoculation loop, touch the surface of a colony of the target species and inoculate the LB
- 3) Incubate at appropriate temperature until $OD_{600}=0.1-0.15$
- 4) Sterilise the glass haemocytometer by spraying with 5% chemogene. Wipe dry and spray again with 70% ethanol. Ethanol should be left to evaporate
- 5) Place a glass cover slip over the centre grid of the haemocytometer
- 6) Pipette 20 μ L of bacterial solution into the trough of the haemocytometer under the coverslip by dispensing at the end of the coverslip
 - a. The liquid should be drawn under the coverslip by capillary action
- 7) Set up brightfield microscope and focus on the haemocytometer
- 8) Using the top left section, count the number of cells in each of the 16 squares
 - a. Do not count cells on the bottom or right hand lines in order to avoid double counting
- 9) Take an average of the 16 squares
- 10) As each small square is 0.1mm³, therefore the number of cells must be multiplied by 10,000 to obtain cells/ml