

# Gas vesicle all lab work

---

MÅNDAG, 2018-06-18

Only Our construct Gas Vesicle 2 has arrived so start with only this one

## **Dilute gblocks and primers**

### **Phusion PCR**

Mixing according to protocol, PCR at 64,8

Result:

Nanodrop DNA concentration measurement : 123,3 ng/μL

### **Nucleic Acid Gel Electrophoresis**

Gel preparation:

- 30mL 1% agarose gel
- 1μL gel red

Gel loading:


1. ladder 1kb
2. sample 1:
  - 1μL dye
  - 0,5μL DNA amplicon
  - 4,5μL MiliQ

Gel run: 95V for 45 min

First try (lane 1 and 2): no result

Second try (lane 4 and 5): still no result

GelDoc-PC\_2018-06-18\_17h01m42s.jpg



Result: failed PCR

Cause: elongation time per cycle was set to 0:30 while a fragment of roughly 2kb requires 2 minutes

ONSDAG, 2018-06-20

#### Digest plasmid (remove MCS)

Cut p413TEF with XbaI and XhoI

- 14µL Nuclease free water
- 2µL 10x fast digest buffer
- 2µL DNA
- 1µL XbaI
- 1µL XhoI

37°C for 5 minutes, no inactivation afterwards

#### Purify plasmid

According to protocol

Nanodrop: 16,6 ng/µL

#### Gradient PCR with Phusion and Dream-Taq


## Phusion PCR preparation

Table1		A	B	C	D	E	F
1	Phusion				Dream-Taq		
2		x1		x9		x1	x9
3	HF buffer		4	36	Dream-Taq buffer	2	18
4	dNTPs	0,4		3,6	dNTPs	0,4	3,6
5	Fw_primer	0,1		0,9	Fw_primer	0,1	0,9
6	Rv_primer	0,1		0,9	Rv_primer	0,1	0,9
7	DNA temp.	0,5		4,5	DNA temp.	0,5	4,5
8	Polymerase	0,2		1,8	Polymerase	0,1	0,9
9	Water	14,7		132,3	Water	18,2	163,8
10	Total		20	180	Total	20	180

## Temperature gradient Phusion and Dream-Taq PCR:

Table2		A	B	C	D	E	F	G	H	I
1			1	2	3	4	5	6	7	8
2	Dream-Taq		60,0	59,4	58,3	56,3	53,9	52,0	50,7	50,0
3	Phusion		70,0	69,5	68,4	66,4	64,0	62,0	60,7	60,0

Result: top: dream-tag, bottom Phusion, 1kb ladder to the left and decreasing gradient from left to right



All Phusion PCRs worked and only 58,3 and 56,3 worked for Dream-Taq

**Purify PCR amplicon**

To maximise concentrations, all 8 Phusion PCR's are collected in 2 tubes: 1-4 in tube 1 and 5-8 in tube 2

Purification according to protocol

nanodrop tube 1: 230,0 ng/μL

nanodrop tube 2: 242,3 ng/μL

**Gibson assembly**P413TEF: 16,6 ng/ $\mu$ LGv2\_1: 230 ng/ $\mu$ L

	A	B
1		Volume
2	p413TEF	5 $\mu$ L
3	Gv2_1	1 $\mu$ L
4	Gibson assembly master mix	10 $\mu$ L
5	MiliQ	4 $\mu$ L
6	Total	20 $\mu$ L

Follow protocol

**E.coli transformation and growth**

Use DH5 alpha cells from -80 freezer

- thaw on ice for 30 min
- add 3 $\mu$ L of assembled plasmid to the cells
- heat shock in 42 $^{\circ}$ C for 45 sec
- put on ice for 20-30 min
- add LB and restore for 1 hour at 37 $^{\circ}$ C
- Streak on plates
- Grow overnight

Result: positive control had growth and negative control not.

Note: negative control was just cells without plasmid, this should have been cells that have transfected with plasmids without construct


TORS DAG, 2018-06-21

**Colony PCR Dream-Taq**

8 colonies were arbitrarily selected from the plate

	A	B
1		Volume (in $\mu\text{L}$ )
2	Dream-Taq buffer	5
3	dNTPs	5 (should have been diluted)
4	Fw_primer	5
5	Rv_primer	5
6	DNA temp.	1
7	Polymerase	1
8	Water	33
9	Total	50

PCR with elongation phase of 2 minutes and annealing temp of 58°C



Result: No construct observed in the PCR

Repeat all the above

Yields plate with transformed E.coli

#### MÅNDAG, 2018-07-02

Incubate two different transformed E.coli colonies at 37°C (t=16:15) in different 10mL tubes

#### TISDAG, 2018-07-03

Take out at 10:00

Do miniprep in duplo on both according to protocol with low copy number plasmid.

Restriction with NcoI and PstI at 37°C for 10 minutes according to Fast Digest Protocol

Gel electrophoresis

Table19	
	A B C D E
1	ladder 1.1 1.2 2.1 2.2

Picture missing

#### TORS DAG, 2018-07-05

Picked 14 colonies of transformed E.coli

Grow both overnight in LB

#### FREDAG, 2018-07-06

Cultures removed from shaker at 10:15

Do mini-prep on all according to protocol

#### Check results:

restriction enzymes NcoI

Gel electrophoresis to check result

No obvious results so restriction and gele electrophoresis is repeated

#### TISDAG, 2018-07-10

Restriction with NcoI and PstI for 10 minutes according to Fast Digest protocol

#### Gel Electrophoresis 2nd try:

Ladder, 3, 4, 5, 6, 9, 10, p413TEF

Run for 40 min at 85V


Gel electrophoresis result: GelDoc-PC 2018-07-10\_14h49m22s

#### FREDAG, 2018-07-13

Restriction with PstI and ScaI for 30 minutes at 37°C according to Fast Digest protocol

loading:

Table17		A	B	C	D	E	F	G	H
1	ladder	1	2	3	4	5	6	uncut ctrl	
2	ladder	8	7	9	10	11	12	cut ctrl	




Need a better separation, repeat with longer run time on long gel tray

MÅNDAG, 2018-07-16

**Gel electrophoresis, longer run**

Table18		A	B	C	D	E	F	G	H
1	ladder	5	6	8	10	11	ctrl cut	ctrl uncut	



Still no good results.

Failed experiment. Repeat from scratch.

TISDAG, 2018-07-24

All primers (listed later) and constructs are diluted to the right concentrations

TORS DAG, 2018-08-02

**Acquire plasmids**

- Get plate with E.coli containing pSP-G2
- Get plate with E.coli containing p414TEF

Do mini-prep on both according to protocol

- Get purified p413TEF
- Get purified p416TEF

Nanodrop DNA concentration measurements:

Table20		A	B	C	D	E	F	G	H	I	J
1	p413	p416	p414_1	p414_2	p414_3	p414_4	PSPG_1	PSPG_2	PSPG_3	PSPG_4	
2	229,8	118,9	27,5	25,4	25,0	6,6	123,4	144,0	120,4	116,5	

TORS DAG, 2018-08-09

**PCR Amplification:**



PTEF/PGK1\_1-3b PCR settings:

- Template DNA: pSP-G2
- Primers: PTEF/PGK1-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 1413 + 2\*15 tail = 1443 bp
- Annealing temperature: 65.3 °C
- Extension time: 1 min 30 sec (approx 1 min per 1kb)

T/ADH1\_1-3b PCR settings:

- Template DNA: pSP-G2
- Primers: T\ADH1\_Fw & -Rv
- Polymerase: Phusion
- Fragment length: 165 + 2\*15 tail = 195 bp
- Annealing temperature: 65.1 °C
- Extension time: 0 min 12 sec (approx 1 min per 1kb, very short, maybe a bit longer to be sure?)

GasVes1 PCR settings:

- Template DNA: gBlock GasVes1
- Primers: Vesic1v2-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 1764 + 2\*15 tail = 1794 bp
- Annealing temperature: 65.5 °C
- Extension time: 1 min 50 sec (approx 1 min per 1kb)

GasVes3b (GvpT) PCR settings:

- Template DNA: gBlock GasVes3b
- Primers: Vesic3b-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 879 + 2\*15 tail = 909 bp
- Annealing temperature: 64,8 °C
- Extension time: 0 min 55 sec (approx 1 min per 1kb)

GasVes2 PCR settings:

- Template DNA: gBlock GasVes2
- Primers: gBlock-Vesci2\_Fw & \_Rv
- Polymerase: Phusion
- Fragment length: 2304 + 2\*15 tail = 2334 bp
- Annealing temperature: 61.1
- Extension time: 2 min 20 sec (approx 1 min per 1kb)

PTEF/PGK1\_KJU PCR settings:

- Template DNA: pSP-G2
- Primers: PTEF/PGK1-GvpU/KJ-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 1413 + 2\*15 tail = 1443 bp
- Annealing temperature: 60.8
- Extension time: 1 min 30 sec (approx 1 min per 1kb)

T/ADH1\_KJU PCR settings:

- Template DNA: pSP-G2
- Primers: T\ADH1-Rv-GvpU & T\ADH1-Fw
- Polymerase: Phusion
- Fragment length: 165 + 2\*15 tail = 195 bp
- Annealing temperature: 65.1

- Extension time: 0 min 12 sec (approx 1 min per 1kb, very short, maybe longer to be sure)

GvpKJ PCR settings:

- Template DNA: gBlock GvpKJ
- Primers: GvpKJ-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 645 + 2\*15 tail = 675 bp
- Annealing temperature: 64.8
- Extension time: 0 min 41 sec (approx 1 min per 1kb)

GvpU PCR settings:

- Template DNA: gBlock GvpU
- Primers: GvpU-Fw & -Rv
- Polymerase: Phusion
- Fragment length: 417 + 2\*15 tail = 447 bp
- Annealing temperature: 64.5
- Extension time: 0 min 25 sec (approx 1 min per 1kb)

First combined PCR at 65,1°C with 01:50 extension time

Do as a gradient PCR between 60-70°C

Table5		A	B	C	D	E	F
1			PTEF/PGK1_1-3 b	<b>GasVes1</b>	<b>GasVes3b</b> (GvpT)	<b>GvpKJ</b>	Overall
2	Template DNA	pSP-G2	<b>gBlock GasVes1</b>	<b>gBlock</b> <b>GasVes3b</b>	<b>gBlock GvpKJ</b>		
3	Primers	PTEF/PGK1-Fw & -Rv	<b>Vesic1v2-Fw &amp;</b> <b>-Rv</b>	<b>Vesic3b-Fw &amp;</b> <b>-Rv</b>	<b>GvpKJ-Fw &amp;</b> <b>-Rv</b>		
4	Annealing T°C	65,30	<b>65,50</b>	<b>64,80</b>	<b>64,80</b>	65,10	
5	Fragment length	1443	<b>1794</b>	<b>909</b>	<b>675</b>		
6	Extension time	01:30	<b>01:50</b>	<b>00:55</b>	<b>00:41</b>	01:50	

Second combined PCR at 64,9°C with 00:25 extension time

Do as a gradient PCR between 60-70°C

	A	B	C	D	E
1		T/ADH_1-3b	T/ADH_KJU	GvpU	Overall
2	Template DNA	pSP-G2	pSP-G2	gBlock GvpU	
3	Primers	T\ADH1_Fw & -Rv	T\ADH1-Rv-GvpU & T\ADH1-Fw	GvpU-Fw & -Rv	
4	Annealing T°C	65,10	65,10	64,50	64,90
5	Fragment length	195	195	447	
6	Extension time	00:12	00:12	00:25	00:25

Third combined PCR at 60,95°C with 02:20 extension time

Do as a gradient PCR between 55-65°C

	A	B	C	D
1		<b>GasVes2</b>	PTEF/PGK1_KJ U	Overall
2	Template DNA	<b>gBlock GasVes2</b>	pSP-G2	
3	Primers	<b>gBlock-Vesic2_Fw &amp; _Rv</b>	PTEF/PGK1-GvpU/KJ-Fw & -Rv	
4	Annealing T°C	<b>61,10</b>	60,80	60,95
5	Fragment length	<b>2334</b>	1443	
6	Extension time	<b>02:20</b>	01:30	02:20

1. Check dilution of PCR primers (10µM)
2. Check dilution of gBlocks (50-250ng needed per tube)
3. Check dilution of plasmid templates (1pg-10ng needed per tube)
4. Fill each PCR tube with:
  - 10µL 5X HF Buffer
  - 1µL dNTPs
  - 2,5µL Fw primer
  - 2,5µL Rv primer
  - 1,5µL DNA template
  - 0,5µL Polymerase
  - fill to 50µL with MiliQ

#### FREDAG, 2018-08-10

Used 5 tubes/samples. Gel results below.


Loading order:

Gel 1, row 1: Ladder, PTEF/PGK1\_1-3b 1-5, GasVes1 1-5, GasVes3b 1-4


Gel 1, row 2: Ladder, GvpKJ 1-5, T/ADH\_1-3b 1-5, T/ADH\_KJU 1-4

Gel 2, row 1: GvpU 1-5, GasVes2 1-5, PTEF/PGK1\_KJU 1-4

Gel 2, row 2: Ladder, GasVes3b 5, T/ADH\_KJU 5, PTEF/PGK1\_KJU 5




📎 20180810\_SuperPcr\_gel2.jpg



**MÅNDAG, 2018-08-13**

Rerun gel for samples GasVes1 3+5, GasVes3b 3-4, GvpKJ 2-5 and GasVes 3-5.

📎 20180813\_reruns\_superpcr.jpg




**TISDAG, 2018-08-14**

Several bands visible for GasVes3b and GvpKJ. 3  $\mu$ l of the third sample of each of these constructs was run on gel and the correct bands were extracted with ThermoFisher Scientific gel extraction kit according to the protocol of the kit.

PCR for GasVes1 and GasVes2 were remade with an extension time of 45 s and an annealing temperature of 64.8 °C for 10 cycles and 72 °C for the remaining 25.

Gel with GasVes1, GasVes2, GasVes3b, GvpKJ.

**ONSDAG, 2018-08-15**

Remade gel extraction of GasVes3b, GvpKJ with 7 \* 3  $\mu$ l of each. Following concentrations were obtained with Nanodrop:

GasVes3b sample 1: 4.3 ng/ $\mu$ l

GasVes3b sample 2: 7.8 ng/ $\mu$ l

GvpKJ sample 1: 14.8 ng/ $\mu$ l

GvpKj sample 2: 9.1 ng/ $\mu$ l

Note that some bands that could be primers were visible on the gel.

PCR purification with ThermoFisher scientific pcr purification kit of samples from 2018-08-09 (PTEF/PGK1\_1-3b sample 4, T/ADH\_1-3b sample 3, T/ADH\_KJU sample 3, GvpU sample 4, PTEF/PGK1\_KJU sample 3 (should double check this)). Final concentrations measured with Nanodrop:

PTEF/PGK1\_1-3b: 69.7 ng/ $\mu$ l

T/ADH\_1-3b: 248.9 ng/ $\mu$ l

T/ADH\_KJU: 199.4 ng/ $\mu$ l

GvpU: 109.9 ng/ $\mu$ l

PTEF/PGK1\_KJU: 73.3 ng/ $\mu$ l



**TORS DAG, 2018-08-16**

Repeated PCR amplification for GasVes1, GasVes2, GasVes3b and GasVesKJ in gradient PCR from 58-68°C

	A	B	C	D	E	F	G	H
1	1	2	3	4	5	6	7	8
2	58,0	58,7	60,0	62,0	64,4	66,4	67,5	68,0

FREDAG, 2018-08-17

Gel electrophoresis of PCR amplicons:



expected bands:

- GasVes1: 1794
- GasVes2: 2334
- GasVes3b: 909
- GasVesKJ: 675


None of the bands match.

new gradient PCR with only GasVes1, GasVes2 and GasVes3b:

Table12									
	A	B	C	D	E	F	G	H	I
		x1	x9		x1	x9		x1	x9
2	HF buffer	4	36	HF buffer	4	36	HF buffer	4	36
3	dNTPs	0,4	3,6	dNTPs	0,4	3,6	dNTPs	0,4	3,6
4	Vesic1v2-Fw	0,1	0,9	Vesic2-Fw	0,1	0,9	Vesic3b-Fw	0,1	0,9
5	Vesic1v2-Rv	0,1	0,9	Vesic2-Rv	0,1	0,9	Vesic3b-Rv	0,1	0,9
6	GasVes1	0,5	4,5	GasVes2	0,5	4,5	GasVes3b	0,5	4,5
7	Phusion Polymerase	0,2	1,8	Phusion Polymerase	0,2	1,8	Phusion Polymerase	0,2	1,8
8	Water	14,7	132,3	Water	14,7	132,3	Water	14,7	132,3
9	Total	20	180	Total	20	180	Total	20	180

2018-08-20

Gel electrophoresis of PCR amplicons:



**Digest plasmid**

Remove MCS: cut plasmids p413TEF, p414TEF and p416TEF

Negative Control: linearise plasmids p413TEF, p414TEF and p416TEF with XhoI

	A	B	C	D	E	F	G
1		p413TEF cut	p414TEF cut	p416TEF cut	p413TEF lin	p414TEF lin	p416TEF lin
2	MiliQ	14	14	14	15	15	15
3	10x FD buffer	2	2	2	2	2	2
4	DNA	2	2	2	2	2	2
5	SacI	0	1	1	0	0	0
6	XhoI	1	1	1	1	1	1
7	XbaI	1	0	0	0	0	0

37°C for 10 minutes (5 min required but 10 to be sure), no inactivation afterwards

### Gel Electrophoresis plasmid digestion:

Gel preparation:

- large tray & 8 well comb
- 30mL 1% agarose gel
- 1µL gel red


Gel sample loading (except ladder):

- 1 µL dye
- 0,5 µL DNA (according to table below)
- 4,5 µL MiliQ

	A	B	C	D	E	F	G	H
1	lane 1	lane 2	lane 3	lane 4	lane 5	lane 6	lane 7	lane 8
2	ladder 1kb	p413TEF cut	p413TEF lin	p414TEF cut	p414TEF lin	p416TEF cut	p416TEF lin	

Gel run: 95V for 45 min

Expected outcome:



### Gel Electrophoresis PCR amplification:

Gel preparation:

- large tray & 16 well comb
- 30mL 1% agarose gel
- 1µL gel red



Gel sample loading (except ladder):

- 1  $\mu\text{L}$  dye
- 0,5  $\mu\text{L}$  DNA (according to table below)
- 4,5  $\mu\text{L}$  MiliQ

Gel run: 95V for 45 min

	A	B	C	D	E	F	G	H	I	J	K
1	lane	lane 1	lane 2	lane 3	lane 4	lane 5	lane 6	lane 7	lane 8	lane 9	lane 10
2	content	ladder 1kb	PTEF/PGK1_1-3 b	T/ADH_1-3b	GasVes1	GasVes2	GasVes3b	PTEF/PGK1_KJ U	T/ADH_KJU	GvpKJ	GvpU
3	expected bands	ladder 1kb	1 at 1443	1 at 195	1 at 1794	1 at 2334	1 at 909	1 at 1443	1 at 195	1 at 675	1 at 447

### DNA Purification Plasmid and PCR amplicons

according to protocol in folder

### Concentration measurement by Nanodrop

	A	B	C	D	E	F	G	H	I	J	K	L	M
1	fragment	p413TEF cut	p414TEF cut	p416TEF cut	PTEF/PGK1_1-3 b	T/ADH_1-3b	GasVes1	GasVes2	GasVes3b	PTEF/PGK1_KJ U	T/ADH_KJU	GvpKJ	GvpU
2	concentration[n g/ $\mu\text{L}$ ]												
3	fragment length	5528	4948	5054	1443	195	1794	2334	909	1443	195	675	447
4	pmols	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!	#VALUE!

ONSDAG, 2018-08-22

### Gibson assembly:

For the tables below, the total Volume of DNA has to be below 10 $\mu\text{L}$ . If not, scale down proportionally to ng quantities to the left.

Assembly of construct GasVes2 into p413TEF:

To calculate:  $m/v = \text{pmols} * \text{bp} * 0,65$

	A	B	C	D
1		pmols	ng	Volume
2	p413TEF cut	0,08	287,46	
3	GasVes2	0,16	242,74	
4	Gibson assembly master mix			10 $\mu\text{L}$
5	MiliQ			10 - V_insert $\mu\text{L}$
6	Total			20 $\mu\text{L}$

Assembly of constructs GvpU and GvpKJ with double promoter and ADH terminator into p414TEF:

	A	B	C	D
1		pmols	ng	Volume
2	p414TEF cut	0,08	257,30	
3	PTEF/PGK1_KJ U	0,16	150,07	
4	T/ADH_KJU	0,42	53,24	
5	GvpKJ	0,16	70,20	
6	GvpU	0,16	46,49	
7	Gibson assembly master mix			10 $\mu$ L
8	MiliQ			10 - V_insert $\mu$ L
9	Total			20 $\mu$ L

Assembly of constructs GasVes1 and GasVes3b with double promoter and ADH terminator into p416TEF:

Table16				
	A	B	C	D
1		pmols	ng	Volume
2	p416TEF cut	0,08	262,8	
3	PTEF/PGK1_1-3 b	0,16	150,1	
4	T/ADH_1-3b	0,42	53,2	
5	GasVes1	0,16	186,6	
6	GasVes3b	0,16	94,5	
7	Gibson assembly master mix			10 $\mu$ L
8	MiliQ			10 - V_insert $\mu$ L
9	Total			20 $\mu$ L

Follow Gibson assembly protocol in folder

Negative control with all 3 cut plasmids and gibson reagents, no constructs.

**Digestion and Gel electrophoresis:** to check succes of Gibson assembly


digest p413 at 37°C with (lane 1): NcoI (in backbone), PstI (in GasVes2)

digest p414 at 37°C with (lane 2): KpnI (in backbone), Bsu15I (in double promoter), XmaJI (in GvpU)

digest p416 at 37°C with (lane 3): KpnI (in backbone), Bsu15I (in double promoter, GasVes1 and GasVes3b)


Expected bands:

- 1: p414 with inserts; Gv2
- 2: p414 with inserts; GvpKJU
- 3: p416 with inserts; Gv1&3b



Gel results:

Very faint bands but seem to indicate succes.



### E.coli transformation and growth

use DH5 alpha cells from -80 freezer

- thaw on ice for 30 min
- add 3µL of Gibson assembled plasmid to the cells
- heat shock in 42°C for 45 sec
- put on ice for 20-30 min
- add LB and restore for 1 hour at 37°C
- Streak on plates
- Grow overnight

**TORS DAG, 2018-08-23**

No growth on the plates, transformation did not succeed

Decision: drop the gas vesicle project to focus on other parts.