

Gel Electrophoresis

Material

Agarose

1xTAE buffer

ddH₂O

6x loading buffer

DNA

GelStain(10000x)

Proceure

1. Prepare 34mL 1xTAE buffer with 1.7g agarose, and boil it three times, then add 1.7 μ L GelStain, shake well.
2. Pour the agarose gel into gel tray, assemble gel pouring apparatus by inserting gate into slots.
3. Allow agarose to cool, place the gel in the apparatus rig with the wells facing the negative end (black-colored).
4. Fill the rig with 1x TAE buffer.
5. Load 2 μ L of DNA maker into lane.
6. Mix 1 μ L of 6x loading buffer with 2 μ L DNA sample, load them into lane.
7. Run at 100V for 30 min.
8. Use the Gel imaging system to check the gel.
9. Take a picture.
10. Throw away the gel carefully and clean up the table-board.