

Paan colour standardization

Protocol:1) To standardize the quantitative assay of colour produced due to mixing of catechu and slaked lime.

Goal: To have a standard graph of absorbance vs concentration of our coloured compound.

The graph obtained shall serve as a standard while doing the enzyme assays to check the efficiency of color degradation.

(KEEPING SLAKE LIME CONSTANT)

Stock: 0.1g% catechu solution

Diluent: Distilled water

Range: 0.01g% - 0.1g%

Absorbance maxima: Needs to be determined.

Ca(OH)₂ :0.1g% aqueous solution

Tube no.	Conc. (g%/ml)	Stock (ml)	Diluent (ml)	Ca(OH) ₂ (ml)	Total Vol.(ml)	Absorbance
1	0.01	0.5	4.5	1.0	6.0	
2	0.02	1.0	4.0	1.0	6.0	
3	0.03	1.5	3.5	1.0	6.0	
4	0.04	2.0	3.0	1.0	6.0	
5	0.05	2.5	2.5	1.0	6.0	
6	0.06	3.0	2.0	1.0	6.0	
7	0.07	3.5	1.5	1.0	6.0	
8	0.08	4.0	1.0	1.0	6.0	
9	0.09	4.5	0.5	1.0	6.0	
10	0.10	5.0	0.0	1.0	6.0	

11	BLANK	----	5.0	1.0	6.0	
----	-------	------	-----	-----	-----	--

In the above table tube no. 2, 4 and 6 will have a catechu: slake lime ratio of 1:1 , 1:2 and 1:3 respectively. (By mass since we are using 0.1g% of the solutions)

(KEEPING CATECHU CONSTANT)

Stock: 0.1g% Slake lime solution

Diluent: Distilled water

Range: 0.01g% - 0.1g%

Absorbance maxima: Needs to be determined. (will be same as above)

Catechu :0.1g% aqueous solution

Tube no.	Conc. (g%/ml)	Stock (ml)	Diluent (ml)	Catechu (ml)	Total Vol.(ml)	Absorbance
1	0.01	0.5	4.5	1.0	6.0	
2	0.02	1.0	4.0	1.0	6.0	
3	0.03	1.5	3.5	1.0	6.0	
4	0.04	2.0	3.0	1.0	6.0	
5	0.05	2.5	2.5	1.0	6.0	
6	0.06	3.0	2.0	1.0	6.0	
7	0.07	3.5	1.5	1.0	6.0	
8	0.08	4.0	1.0	1.0	6.0	
9	0.09	4.5	0.5	1.0	6.0	
10	0.10	5.0	0.0	1.0	6.0	
11	BLANK	---	5.0	1.0	6.0	

Protocol:2) To measure the optical density of the colour produced due to mixing of catechu and slaked lime colorimetrically.

(KEEPING CATECHU CONSTANT)

- 1) Weigh fixed amount of sieved catechu powder and slaked lime (in grams) (amount needs to be standardised)
- 2) Mix the two components in various ratios.(1:2,1:3,1:4,1:5) (since catechu is constant the amount of slaked lime must be varied)
- 3) Add a fixed volume of boiling water to the mixtures with 1:2,1:3,1:4 and 1:5.

Ratio	Amt of catechu (in g)	Amt of slaked lime (in g)	Boiling water (in ml)	Absorbance
1:1	0.1	0.1	100	
1:2	0.1	0.2	100	
1:3	0.1	0.3	100	
1:4	0.1	0.4	100	
1:5	0.1	0.5	100	

(KEEPING SLAKE LIME CONSTANT)

- 1) Weigh fixed amount of sieved catechu powder and slaked lime (in grams) (amount needs to be standardised)
- 2) Mix the two components in various ratios.(1:2,1:3,1:4,1:5) (since slaked lime is constant the amount of catechu must be varied)
- 3) Add a fixed volume of water to the mixtures with 1:2,1:3,1:4 and 1:5.

Ratio	Amt of slaked lime (in g)	Amt of catechu (in g)	Boiling water (in ml)	Absorbance
1:1	0.1	0.1	100	
1:2	0.1	0.2	100	

1:3	0.1	0.3	100	
1:4	0.1	0.4	100	
1:5	0.1	0.5	100	

Protocol:3) To determine the relative amounts of Catechu and slake lime required in media preparation.

For Catechu+slaked lime:

- 1) Take six 100 ml beakers and add 100 ml of boiling distilled water to it.
- 2) To each of the beakers add 0.1g gram of catechu and let it boil again for 10 minutes to extract the catechu pigment.
- 3) Cool these solutions and then add slake lime paste to it in the following manner:

Ratio (Catechu:slake lime)	Amount of catechu(grams)	Boiling water (ml)		Slake lime (grams)	Absorbance
1:0.1	0.1	100		0.01	
1:0.2	0.1	100		0.02	
1:0.3	0.1	100		0.03	
1:0.4	0.1	100		0.04	
1:0.5	0.1	100	Cool the	0.05	
1:0	0.1	100	solution	0.0	
1:1	0.1	100	For 15-20	0.1	
1:1.5	0.1	100	mins	0.15	
1:2	0.1	100		0.2	

Protocol:4) To determine the effect of varying concentration of catechu extracted and constant amount of slake lime added.

From the experimentation done with protocol 4 it was determined that maximum intensity of red colour solution was obtained in solution of 1:0.1 catechu:slake lime concentration.

- 1) Boil around 500 ml of distilled water and using a measuring flask transfer 100 ml to 5 different beakers.
- 2) Add 0.1, 0.2, 0.3, 0.4, 0.5 grams of fine catechu powder to different beakers.
- 3) Hence ratios would now be, 10:1, 20:1, 30:1, 40:1, 50:1 catechu:slake lime conc.
- 4) Heat these mixtures again separately till boiling and let cool at RT then ice water bath.
- 5) Add 0.01 g slake lime in each beaker and then stir till the white powder is completely dissolved.
- 6) Centrifuge the mixtures and collect the supernatant.
- 7) Discard the pellet (which may be white ppt or slightly reddish brown)

Protocol: 5) To determine the relative amounts of Catechu and slake lime required in media preparation.

Goal: To have a standard protocol for extraction of catechu pigment for the purpose of media preparation.

- 1) Take six 100 ml beakers and add 100 ml of boiling distilled water to it.
- 2) To each of the beakers add 0.1g gram of catechu and let it boil again for 10 minutes to extract the catechu pigment.
- 3) Cool these solutions and then add slake lime paste to it in the following manner:

Ratio (Catechu:slake lime)	Amount of catechu(grams)	Boiling water (ml)		Slake lime (grams)	Absorbance
1:0.1	0.1	100		0.01	
1:0.2	0.1	100		0.02	
1:0.3	0.1	100		0.03	
1:0.4	0.1	100		0.04	
1:0.5	0.1	100	Cool the	0.05	
1:0	0.1	100	solution	0.0	
1:1	0.1	100	For 15-20	0.1	
1:1.5	0.1	100	mins	0.15	
1:2	0.1	100		0.2	