

IMPACT-CN

- Preparation of Plasmids
 - Clone your protein fragment into one of the following vectors:
 - pTXB1
 - pTYB1
 - pTYB2
 - pTYB11
 - pTYB12
 - transfer the complete construct into ER2566 via electro transformation
- cell growth
 - prepare a pre-culture with some cell material within up to 10mL LB substituted with 200 µg/mL Amp and 20 µg/mL Cm
 - incubate the cultures overnight at 37°C
 - prepare the main culture:
 - 250mL LB (+Amp and Cm)
 - Inoculate the culture to an OD₆₀₀ of 0.1 with the pre culture
 - Incubate the main culture at 37°C until an OD₆₀₀ auf 0.5 to 0.8 is reached
 - Induce the protein expression via addition of IPTG (0.3-0.5mM)
 - (1M stem solution: 125 µL per culture)
 - Let the cells grow further for 30min at 37°C
 - Protein expression over night at 19°C
- Cell harvest and decomposition

- Pellet the whole culture in large centrifugal jars
 - 4,500rpm, 10-20min, 4°C
- Discard the supernatant
- Resuspend the pelleted cells in 20-30mL lysis solution
 - Retransfer them into falcon tubes
- Decompose the cells using the French Press
 - 2 times for each sample
 - Pressure of approximately 1200
 - Work with ice whenever possible
- Centrifuge the samples to divide the desired substances from the cell debris
 - 4,500rpm, 1h, 4°C
 - Transfer the supernatant into falcon tubes
- Column preparation
 - Close the lower part of the columns with Para film
 - Fill them with 5mL of water(=1 volume) and mark the fill level
 - Let the water out
 - Fill the column with chitin beads up to the marking
 - An evenly surface gets better results
 - Put a filter plate onto the beads
 - **From here, all steps should be carried out at 4°C**
 - Equilibrate the columns with 10volumes of column buffer (low-salt)
- Loading of the columns
 - Put the supernatant onto the columns
 - Make sure the flow-through is less than 0.5-1 mL/min
 - E.g. with the use of clamps
- Washing of the columns
 - Wash the columns with 10 volumes of the column buffer (high salt)
- Induction of cleavage

- Wash the columns with 3 volumes of cleavage buffer
 - Fast, so no clamps!
- Let the columns induce at 4-23°C for at least 16 (up to 40) hours
- Elution
 - Eluate the columns with 1-2 volumes column buffer (low salt)
 - Collect the eluate within falcon tubes or blue caps
 - Concentrate the eluate by using suitable protein filters
 - Centrifuge for 30-45 minutes at full speed up to an end volume of about 250 μL
 - Wash the concentrated protein solution with 4mL protein wash buffer
 - Repeat this 2-3 times total
 - Optionally: wash with storage buffer to ensure a longer functionality
 - Make an Bradford Estimation
 - E.g. with Roti® Nanoquant(Roth), see our other protocols
 - Store the proteins at the freezer or refrigerator
 - Depending on the usage
- Column regeneration
 - All of the columns can be reused 4-5 times
 - Resolve residual proteins with 3 volumes of NaOH (0.3M)

- Incubate for 30minutes
- Wash the columns with further 7 volumes NaOH
- Wash with 20 volumes of H₂O
- Wash with 5 volumes of column buffer
- Store the columns at 4°C, closed with Para film

- Solutions needed for this protocol:

- Lysis Buffer

Substance	concentration	mass/volume in 1L
Triton X-100	0.1%	1mL
PMSF	20 µM	200 µL
TCEP	0.1-1mM	100-1000 µL
Make the following stem solutions: <ul style="list-style-type: none"> • 100mM PMSF • 1M TCEP Resolve the substances within column buffer (low salt)		

- Column Buffers

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3,6g
	<ul style="list-style-type: none"> • 2 H₂O 	
NaCl	50mM (low salt)	2,92g
	500mM (high salt)	29,22g
	1000mM(high salt)	58,44g
	1500mM(high salt)	87,66g

EDTA	1mM	10mL of stem solution
<ul style="list-style-type: none"> • Prepare an EDTA stem solution of 0.1M, resolve in ddH₂O • Resolve all substances in ddH₂O • Adjust the pH to 8.0 (6.0-9.0) with phosphoric acid 		

○ Cleavage Buffer

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3.6g
NaCl	50mM (optionally up to 1000mM)	29g
EDTA	1mM	10mL
DTT	50mM	*
<ul style="list-style-type: none"> • *0.93g DTT resolved in 6mL ddH₂O+119mL column buffer (low salt) make up 125mL cleavage buffer • All substances resolved in H₂O • pH has to be adjusted to 7-9 with phosphoric acid • Store dark because of the DTT! 		

○ Protein wash buffer

substance	concentration	mass/volume in 1L
Na ₂ HPO ₄	20mM	3.6g
NaCl	10mM	0.58g
EDTA	1mM	10mL

From: iGEM Bielefeld-CeBiTec