

LacI-P_{lac}

●Maximize the repression effect

1. Purpose of the experiment:Maximize the repression effect of LacI protein.

1.1. Reagent: The objective strain(pGLO-LacI) , Amp, liquid minimal medium,arabinose,fructose.

1.2. Instrument: liquid transfer guns and gun heads,tubes (each containing 5mL liquid minimal medium), EP tube, 37 °C oscillating incubator, etc.

2. Experimental steps

①Strain contained target plasmid was cultured in the test tubes with liquid minimal medium. And

the concentration gradients of arabinose were set in these tubes (0M, 2×10^{-5} M, 8×10^{-5} M) .

Number	20mM fructose/ μl	20mM arabinose/μl	Amp/μl
0	80	0	5
1	80	10	5
2	80	40	5

②OD₆₀₀ and GFP values of the bacterial culture were measured periodically with the automatic microplate reader. For green fluorescence, excitation and emission wavelengths were 395 and 509 nm.

③Draw charts using measured data to observe the trends of repression.

3. Experimental result




Fig.1 The Relative Fluorescence Unit of the objective strain

4. Details

[illegible]

8.1 10: 00 (23h)										
classification	ara(10mM) /μL	OD ₆₀₀			Fluorescence when OD ₆₀₀ is 0.4			Relative Fluorescence Unit		
		1	2	Average	1	2	Average			
pGLO	0	0.525	0.557	0.541	781.37	787.10	784.24	269.32		
	10	0.432	0.451	0.442	1856.90	1912.80	1884.85	1369.93		
	40	0.472	0.475	0.474	1997.50	2058.10	2027.80	1512.88		
pGLO+lacI	0	0.554	0.573	0.564	1001.10	921.73	961.42	446.50		
	10	0.474	0.493	0.484	681.36	720.71	701.04	186.12	260.39	0.3962
	40	0.399	0.411	0.405	686.13	703.14	694.64	179.72	266.79	0.4059
PUC18	0	0.498	0.528	0.513	508.97	520.87	514.92	0.00		
	10	0.437	0.450	0.444	498.10	538.13	518.12	/		

[illegible]

[illegible]