

(1) Inoculated culture

1. The purpose of the experiment: prepare for the extraction of plasmids

(2) Extraction plasmid

1. The purpose of the experiment: prepare for PCR

2. Experimental steps:

5ml overnight cultured liquid was taken and centrifuged at  $8000 \times g$  for 2 minutes to collect the bacteria and discard the culture medium. Adding 250 $\mu$ l Buffer P1 to the precipitate, suspending the bacteria thoroughly. Add 250 $\mu$ l Buffer P2, immediately gently reverse EP tube 5-10 times, leave at room temperature for 2-4 minutes. After adding 350 $\mu$ l Buffer P3, the EP tube was reversed for 5-10 times and centrifuged at  $12000 \times g$  for 5-10 minutes. Transferred the supernatant fluid to the adsorbent column and centrifuged at  $8000 \times g$  for 30 seconds, and the liquid in the collecting tube was emptied. Add 500 $\mu$ l Buffer DW1, centrifuged at  $9000 \times g$  for 30 seconds and the liquid in the collecting tube was emptied. Add 500 $\mu$ l Wash Solution, centrifuged at  $9000 \times g$  for 30 seconds and the liquid in the collecting tube was emptied. Add 500 $\mu$ l Wash Solution, centrifuged at  $9000 \times g$  for 30 seconds and the liquid in the collecting tube was emptied. The empty adsorption column was centrifuged for 1 minute at  $9000 \times g$ . The adsorption column was placed in a clean 1.5ml EP tube and 30 $\mu$ l ddH<sub>2</sub>O was added in the center of the adsorption membrane. Centrifuged at  $9000 \times g$  for 1 minute. Then the liquid in EP tube was sucked out by liquid transfer gun, and added to the center of adsorption membrane and added 30 $\mu$ l ddH<sub>2</sub>O. After leaving at room temperature for 1 minute, centrifuged at  $9000 \times g$  for 1 minute, and the plasmid DNA solution in the tube was preserved.

PS: 1. ddH<sub>2</sub>O preheated in the oven at 55 °C.

(3) measurement of plasmid sample concentration

1. The purpose of the experiment: to provide the basis for determining the formula of PCR system.

2. experimental steps:

1 $\mu$ l ddH<sub>2</sub>O was added to the measuring platform with a liquid transfer gun, and the DNA concentration was measured (when the concentration was  $> 0.1 \text{ ng}/\mu\text{l}$ , the measurement platform was not cleaned and ddH<sub>2</sub>O was re-measured until the concentration was lower than 0.1  $\text{ng}/\mu\text{l}$ ). The DNA concentration of 1 $\mu$ l solution was measured and the concentration was recorded.

3. experimental result

	dusk-Cas9-pUC57①
DNA concentration ( $\text{ng}/\mu\text{l}$ )	58.2692