



New England BioLabs-Transformation protocol for DH5 α

- 1) Thaw 50 μ l of DH5 α chemically competent cells on ice.
- 2) Mix cells with 1 μ l of template. Gently pipette up and down a few times. Incubate tube on ice for 30 minutes.
- 3) Heat shock tubes at 42°C for 30 seconds.
- 4) Place on ice for 5 minutes.
- 5) Add 950 μ l of SOC.
- 6) Incubate at 37°C for 1 hour shaking at 200 rpm.
- 7) Plate 100 μ l on a petri plate with LB media and 25 mg/ml chloramphenicol.
- 8) Spin down cells at full speed for 3 minutes. Discard 800 μ l of the supernatant.
Resuspend the cells in the remaining 100 μ l. Plate the 100 μ l on a petri plate with LB media and 25 mg/ml chloramphenicol.