Freezing and Thawing of SHSY5Y

How to freeze SHSY5Y for storage at -80°C and reconstitute them.

Materials

- Basic Culture Media
- DMSO
- Centrifuge Tube
- Freezing vials

Procedure

Table 1 Freezing Media (18 ml)

Concentration	Chemicals
17,1 ml	Basic Culture Medium
900 μ1	DMSO

Medium should contain 5% (v/v) DMSO

Split cells and count. One freezing veil should contain somewhere around 2-6*10^6 cells. Calculate the number of cells and quantity of desired vials.

- Pre-cool veils and medium to 4°C
- Transfer required amount of cells to a centrifuge tube
- Spin down for 2 minutes at 1.000 x g and 4°C.
- Discard supernatant
- Resuspend cells in 1ml Medium per desired vial
- Transfer 1 ml resuspend cells in each vial
- Keep vials on 4°C for 10 minutes, then at -20°C for further 10 minutes
- Keep cells at -80°C for 24 hours
- Transfer for final storage to liquid nitrogen

Thawing

Equipment:

- Basic Culture Medium
- Vial of frozen SHSY5Y
- Centrifuge tube
- 75 cm² culture flask
- Take vial from liquid nitrogen and keep it on ice for 10 minutes
- Rapidly thaw in 37°C water bath until defrosted
- Transfer to 9 ml room temperature Basic Culture Medium
- Invert gently
- Centrifuge for 2 minutes at 1.000 x g
- Discard supernatant
- Resuspend in 10 ml Basic Culture Medium
- Transfer to 75cm² flask
- Grow the cells in the Basic Culture Medium at 5% CO2 and 37°C.
- Change medium every 2-3 days for maintaining SHSY5Y cells in culture

References

Shipley, Mackenzie M., Colleen A. Mangold, and Moriah L. Szpara. "Differentiation of the SH-SY5Y human neuroblastoma cell line." *Journal of visualized experiments: JoVE* 108 (2016): 53193.