

This protocol is derived from QIAGEN's RNeasy Protect Bacteria Mini Kit [1].

## Experiment

### Materials

- QIAGEN RNeasy Mini Kit Protect Bacteria - containing:
    - RLT buffer
    - TW1 buffer
    - RPE buffer
    - Purification columns
    - RNase-free collection tubes.
    - RNAProtect Bacteria Reagent
  - Two cultures of *E.coli*, approx. 6ml each, grown in M9-media to an average  $OD_{600}$  of 0,75. One culture is grown alongside with live strongyles while the other one isn't.
  - TE-buffer (1M Tris-HCl, 0,1M EDTA)
  - Lysozyme (s)
  - Proteinase K (20mg/ml)
  - 2-mercaptoethanol (99%)
  - Ethanol (99%)
  - RNase-free water
  - DNase enzyme (1500U)
  - Centrifuge
  - Pipettes, tips and appropriate tubes
- 
- All relevant buffers and reagents were prepared to working solutions:
  - TE-buffer - 2ml, mixed with 15mg/ml lysozyme, 20µl/ml Proteinase K
  - RLT-buffer - 2,5ml, mixed with 10µl/ml 2-mercaptoethanol
  - RPE-buffer - 11ml, four volumes of ethanol (99%) are added to a total of 55ml.
  - DNase (diluted with RNase-free water to 3U/µl)

### Method

1. Three samples are taken from each *E.coli* culture - one aliquot of 1,6ml, 1,7ml and 1,8ml. This is to measure any difference in RNA yield or quality when performing the experiment with different cell amounts.

- 2.** RNeasy Protect Bacteria Reagent is added to each aliquot to a volume totaling 3x the original volume.
- 3.** The aliquots are centrifuged for 10min at 13.000rpm. Supernatant discarded, only pellet remains.
- 4.** 200µl of TE-buffer were added to the samples.
- 5.** Samples were vortexed briefly and incubated for 5 minutes with shaking at 25 degrees celsius.
- 6.** 700µl of RLT-buffer were added to each sample. The samples were centrifuged at 13.000rpm for 2min to remove particulate matter - the supernatants were moved to new tubes.
- 7.** 500µl of ethanol were added to each sample.
- 8.** 700µl of the samples were moved over to RNeasy spin columns (six columns in total).
- 9.** 350µl of RW1 buffer were added to each sample. The columns were centrifuged for 15sec at 13.000rpm. Flowthrough discarded.
- 10.** 10µl of DNase solution were mixed with 70µl of RDD-buffer. 80µl of the DNase/RDD buffer was added to each sample. The samples were incubated for 15min at 25c.
- 11.** 700µl of RW1 buffer were added to each sample. Incubated for 5min at 25c. The samples were then centrifuged for 15sec at 13.000rpm. Flowthrough and collection tubes discarded; new collection tubes attached to the spin columns.
- 12.** 500µl of RPE-buffer were added to each sample. The samples were centrifuged for 15sec at 13.000rpm. Flowthrough discarded.
- 13.** 500µl of RPE-buffer were added to each sample. The samples were centrifuged for 2min at 13.000rpm. Flowthrough and collection tubes discarded. New collection tubes were attached to the spin columns.
- 14.** 30µl of RNase-free water were added to each column. The columns were centrifuged for 1min at 13.000rpm. The flowthrough was collected and put on ice - it contains purified RNA.
- 15.** 15µl of RNase-free water were added to each column to collect any remaining RNA. The columns were centrifuged for 1min at 13.000rpm. The flowthrough was collected and put on ice.

**16.** The RNA suspensions (twelve in total) were measured with Thermo Fishers NanoDrop to check the nucleic acid concentration and sample purity. A gel electrophoresis was also performed using the samples collected in step 15.

## **References**

[1] QIAGEN, 2018. RNAProtect Bacteria Reagent Handbook.

<https://www.qiagen.com/us/resources/resourcedetail?id=95346297-ae79-41a8-9259-e48f702d4f36&lang=en> Date of Visit 2018-10-15