

Protocol Name: Re-isolating endophytes.

Category: Endophyte development

Date: 03/10/2018

Author: Lewis Tomlinson

Source(s): Arnold, A. E. & Lutzoni, F. (2007) Diversity and host range of foliar endophytes: Are tropical leaves biodiversity hotspots? Ecology 88, 541-549

Time Required: Approx. 10 mins

Additional Notes: This protocol is adapted from Arnold and Lutzoni (2007) for use with young Arabidopsis seedlings.

Materials

- Agar plates relevant to the desired endophyte (TSA in our case).
 - Plants (Arabidopsis seedlings in our case)
 - Spreaders
 - NaOH (50%)
 - Ethanol (95-100%)
 - Ethanol (70%)
 - Petri dishes
 - water
-

Procedure.

- Gently rinse plant with water.
- Place plant in a petri dish containing 95% ethanol for 10 seconds and gently shake.
- Place plant in petri dish containing 50% NaOH and gently shake for 1 minute.
- Place in petri dish containing 70% ethanol for 1 minute.
- In a sterile environment allow plant to slightly dry before cutting into desired segments.
- Using spreaders gently smear the extrusion from the plant tissue onto the growth media and incubate at the optimum temperature of your desired endophyte.