

## 06. (June) 2019

---

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Theresa Keil

THURSDAY, 13/6/2019

---

centrifuge exosome supernatant 20 min at 4,000 x g to pellet cells

His-Tag purification using protocol form 06.06.19, adjusted pH of elution buffer to pH 8.0

labeld samples S - supernatant 13.06.19

w - wash

e1 - elution 1

e2 - elution 2

e3 - elution 3

recycling of Ni-NTA resing according to His-TAG-Purification protocol

IVT-fluc-MS2 gel electrophoresis

bands were in correct sizes (1.8 kb)

Gel purification with gel dissolving buffer (DNA Gel CleanUP Kit) - 4x the initial volume

continue on step 4 from RNeasy mini Kit (transfer 700 µl sample to the column)

After dissolving buffer we also added DNase

Concentration 2.0 ng/µl