



Lethbridge HS iGEM 2019

Size Exclusion Chromatography Purification of Cas13a

Buffers

S200 Buffer:

10 mM HEPES (pH 7.5)

1 M NaCl

5 mM MgCl_2

2 mM DTT

Storage Buffer:

600 mM NaCl

50 mM Tris-HCl (pH 7.5)

5% v/v glycerol

2 mM DTT

Column Preparation

Purified using a Superdex75 large column from GE

- 1) Hook column up to AktaPrime chromatography system.
- 2) Wash column with at least 3 column volumes of MilliQ water (900 mL) at a rate of 1.5-2.0 mL per min.
- 3) Equilibrate column in 3 column volumes of S200 buffer (900 mL) at a rate of 1.5 mL per min.

Purification Procedure

- 1) Wash the loop with 10 mL of MilliQ water, followed by 10 mL of S200 buffer.
- 2) With system set to load, inject sample into the loop.



Lethbridge HS iGEM 2019

- 3) Set manual program to monitor UV at 215 nm, 254 nm, and 280 nm, and a 1 MPa pressure alarm. Switch to inject sample onto the column, and fractionate 2 mL fractions at a rate of 1.0 mL per min. Execute.
- 4) Collect fractions for about 250 mL, save fractions that have an A_{280} peak. Take a 50 μ L sample for SDS PAGE analysis and store at 4 °C.
- 5) If another sample is being loaded onto the column, wash loop with 10 mL of buffer first before loading sample.
- 6) After purification(s) have finished, wash column with two column volumes (600 mL) of MilliQ water, followed by two column volumes of 20% ethanol (600 mL) for storage.

Concentration

- 1) Rinse Vivaspin MWCO 30000 with storage buffer (2 mL) and centrifuge at 4000 xg for 10 min
- 2) Pool elutions together.
- 3) Remove buffer from Vivaspin and add pooled elutions.
- 4) Centrifuge at 4000 xg for 5 min
- 5) Keep centrifuging: check for speed of concentrating. Concentrate to 1-2 mL. Watch for precipitate!
- 6) Remove filtrate from bottom of Vivaspin, and pool with other filtrate and store at 4 °C.
- 7) Aliquot protein sample and flash freeze in liquid nitrogen and store at -80 °C until needed.