

# CHEMOCOMPETENT CELLS PREPARATION AND TRANSFORMATION

ADAPTED FROM THE BONNET TEAM PROTOCOL REPOSITORY

## MATERIALS:

### To prepare 50mL of culture :

- LB medium
- Antibiotics (according to the plasmid used for transformation)
- TSS medium: To make 50 mL : 5g PEG 8000, 1.5 mL 1M MgCl<sub>2</sub> (or 0.30g MgCl<sub>2</sub>\*6H<sub>2</sub>O), 2.5 mL DMSO and LB to 50 mL. Filter sterilized (0.22 µm filter).
- LB plates
- 1.5mL tubes
- 50mL falcons
- Ice
- 500mL flask for culture
- Liquid nitrogen

### For transformation:

- Chemical competent cells
- Ice
- Clean DNA to transform
- SOC
- Selective agar plates

## PROTOCOL:

### Cells preparation

Day 1:

- Dilute the overnight culture into 50mL of LB without antibiotic at 1/500 (200µL in 100mL)

- Incubate at 37°C until OD<sub>600</sub> reaches 0.2-0.3 (3-4h)
- Place 1.5mL tubes, racks, 10mL pipettes at -20°C
- **FROM NOW ON DO EVERYTHING ON ICE**
- Incubate the culture on ice for 10min in 50mL falcon tubes
- Cool down the centrifuge to 4°C
- Centrifuge the culture at 3000rpm 4°C for 10min
- Remove the supernatant
- Resuspend cells in 10% volume of TSS buffer
- Aliquot cells in 100uL in 1.5mL tubes with Multipette
- Flash freeze them in liquid nitrogen
- Store at -80°C

**To test them: use pUC19 as positive control and do not forget negative control.**

- Thaw gently 3 tubes of cells on ice
- Add 1μL pUC19 in 1 tube, add nothing in the others (negative control) (keep cells on ice)
- Incubate 30 min on ice
- Heat-shock cells at 42°C during 45s (in water bath)
- Put back on ice for 2 minutes after heat-shock
- Add 900 μL pre-warmed SOC (37°C) (rich medium)
- Incubate cells at 37°C with agitation during at least 30min
- Centrifuge cells at 4000rpm for 2min, remove 800μL of supernatant
- Plate the rest from pUC19 positive control on LB Carb plate
- Plate 100μL from negative control respectively on LB Chloramphenicol, LB Kanamycin, LB Carbenicillin and LB Spectinomycin plates
- Incubate at 37°C overnight

Efficiency should be around  $10^7$  colonies/ $\mu\text{g}$  pUC19

**Transformation protocol:**

- Thaw gently cells on ice
- Add  $1\mu\text{L}$  DNA (keep cells on ice)
- Incubate on ice for 30 min
- Heat-shock cells at  $42^\circ\text{C}$  during 45s (in water bath)
- Put back on ice after heat-shock for 2 min
- Add  $900\mu\text{L}$  pre-warmed SOC ( $37^\circ\text{C}$ ) (rich medium)
- Incubate cells at  $37^\circ\text{C}$  with agitation during at least 30min
- Plate  $100\mu\text{L}$  of transformation on selective LB agar plate or centrifuge cells at 4000rpm during 1min, remove  $800\mu\text{L}$  of supernatant and plate the resuspended pellet
- Incubate at  $37^\circ\text{C}$  overnight